zimmer

31 luglio 2013

A: Chirurghi

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22

Trabecular Metal™ Reverse Shoulder System Surgical Technique

Implant Insertion

Glenosphere Assembly

The Glenosphere is typically inserted prior to humeral component final seating to maximize exposure of the glenoid and ease of insertion. Ensure all soft tissue is removed around the Base Plate to allow the Glenosphere to completely seat.

Assemble the Glenosphere Helmet Inserter by threading the Dual Taper/Spacer Impactor into either the 36mm (green) or the 40mm (yellow) Glenosphere Helmet (Fig. 71). Insert the appropriate diameter Glenosphere into the helmet by sliding it into the helmet so that the Glenosphere is held in place by the body of the helmet and the tabs rest securely underneath the Glenosphere (Figs. 72 & 73). Wipe the Base Plate taper clean of all fluids and inspect the taper to ensure it is free of scratches or damage. Place the Zimmer Shoulder Shoehorn Retractor on the posterior side of the glenoid to aid in retracting the humerus and other soft tissue (Fig. 72). When approaching the Base Plate, a finger can be placed on top of the Glenosphere to help guide and feel the Glenosphere slide over the taper into position.

Note: While engaging the Glenosphere, it is important to monitor the position of the proximal humerus and provisional along with retractors since they could interfere with Glenosphere placement. Alternatively, a bone hook can be placed on the humeral provisional to draw the humerus laterally to provide clearance for the glenosphere. It's important to feel the mechanical resistance of taper engagement before proceeding to impaction.

Once the Glenosphere is seated evenly and circumferentially, use your free hand to press firmly on the Glenosphere to secure it to the Base Plate. Keeping

a finger on the Glenosphere, remove the Glenosphere Helmet pulling the instrument away in the SAME DIRECTION used to insert the Glenosphere (i.e. If an anterior approach was used to insert the Glenosphere, remove the instrument by pulling it from the anterior direction).

This will help minimize changes to the Glenosphere placement on the Base Plate and damage to the Glenosphere Helmet itself.

Fig. 71


Fig. 72


Fig. 73


Fig. 74

Note: If unable to visually confirm an even, circumferential engagement of the Glenosphere to the Base Plate, consider the use of a fluoroscope to aid in the confirmation. Seating of the Glenosphere to the Base Plate can be examined in the axillary view or in a view parallel to glenoid version. The medial rim of the glenosphere should be parallel to the face of the Base Plate (Fig. 75).

Assemble Glenosphere Impactor Head to the Impactor Handle and place the Glenosphere Impactor Head centrally on the Glenosphere. Strike the Glenosphere Impactor Head with 3 firm mallet strikes to engage the Glenosphere on the Base Plate (Fig. 76). Pull on the Glenosphere to verify the taper is locked. Reconfirm uniform engagement between the Base Plate and Glenosphere by using a small angled or 90 degree clamp to assess for malalignment gaps anterior to posterior, as well as inferior to superior.


Fig. 75


Fig. 76

Glenosphere Removal

Should it become necessary to remove the Glenosphere, the Glenosphere Distractor can be used. Assemble the Glenosphere Distractor. Wedge the fin tip between the superior glenoid bone and the underside of the Glenosphere (Fig. 77). There must be good contact on these two surfaces for disengagement to occur. Pull the Glenosphere Distractor trigger until it fires. The Glenosphere head should be loose enough to gently remove by hand. If not, repeat the step making sure there is contact between the distractor tip, the glenoid bone surface and the Glenosphere head. Trial if necessary and implant the final Glenosphere as described on pages 19-21. Reduce the joint, and confirm range of motion. If all is satisfactory, continue on to the Closure Section.

Fig. 77

Humeral Stem Insertion

The final humeral preparation for a 130mm length Trabecular Metal humeral stem may be identical to a 130mm length Non-Porous Reverse humeral stem, allowing intra-operative flexibility between the two as a final implant. This includes preparation for the 10, 12, 14, 16 and 18mm sizes. When implanting an 8mm humeral stem, there are distinct proximal reamers for the two respective stems. The 8mm Non-Porous Reverse proximal reamer is marked with an etch reading "8MM NP" on the reamer shaft to differentiate it from the proximal reamer used for an 8mm Trabecular Metal humeral stem. The proximal geometry of the 8mm Non-Porous Reverse humeral stem is smaller than the 8mm Trabecular Metal humeral stem.

Cemented Technique Humeral Preparation

If using a Cement Restrictor Plug, insert the plug one centimeter distal to the tip of the Humeral Stem. Thoroughly clean and dry the canal. Inject cement into the humeral canal. Use a finger to thoroughly pack the cement.

Note: Stem size is chosen based on cement mantle desired and the last reamers used.

Technique Tip: Be careful to avoid contact between the cement and the Trabecular Metal material as the cement will interfere with the biological ingrowth properties of the material.

Note: To avoid risk of a periprosthetic fracture, ensure the final humeral stem implant is not larger than the last Intermedullary Reamer used. Similarly, to avoid an overly loose canal fit, ensure the final humeral stem implant is not smaller than the last Intramedullary Reamer used.

Press-fit Technique

The Humeral Stem can be press-fit by sizing to the reamed diameter. Refer to sizing chart on page 4 for press-fit/clearance conditions.

Insertion for Cemented and Press-fit Techniques

Before inserting the final humeral component, drill any desired suture holes through the proximal neck of the humerus. Attach the Humeral Stem Inserter/Extractor to the Humeral Stem Implant by opening the handle all the way, inserting the stem inserter end into the proximal opening in the assembly, and closing the handle to lock the inserter in place (Fig. 78).


Fig. 78