25 settembre 2014

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INTRO	Surgical Technique	
<div style="margin-bottom: 20px;"> <h3 style="color: #0070C0;">Magnet Usage</h3> <p>Warning: Some instruments in the <i>Persona</i> System contain magnets. All <i>Persona</i> magnetic instruments should be kept at a safe distance from a patient's active implantable medical device(s) (i.e. pacemaker). These types of devices may be adversely affected by magnets. Instruments containing magnets should be kept on an appropriate table or stand when not in use at the surgical site.</p> </div> <div style="margin-bottom: 20px;"> <h3 style="color: #0070C0;">Symbols</h3> <p>Symbols have been established for the following:</p> <ul style="list-style-type: none"> Left Right Varus/Valgus Medial/Lateral Standard Do not implant – Not for implant Lock Unlock Anterior Referencing Do not impact Cemented Stemmed Inset Only </div> <div style="display: flex; flex-wrap: wrap; justify-content: space-around;"> <div style="text-align: center; margin: 5px;">  <p>Left</p> </div> <div style="text-align: center; margin: 5px;">  <p>Right</p> </div> <div style="text-align: center; margin: 5px;">  <p>Varus/Valgus</p> </div> <div style="text-align: center; margin: 5px;">  <p>Medial/Lateral</p> </div> <div style="text-align: center; margin: 5px;">  <p>Standard</p> </div> <div style="text-align: center; margin: 5px;">  <p>Do not implant - Not for implant</p> </div> <div style="text-align: center; margin: 5px;">  <p>Lock</p> </div> <div style="text-align: center; margin: 5px;">  <p>Unlock</p> </div> <div style="text-align: center; margin: 5px;">  <p>Anterior Referencing</p> </div> <div style="text-align: center; margin: 5px;">  <p>Do not impact</p> </div> <div style="text-align: center; margin: 5px;">  <p>Cemented</p> </div> <div style="text-align: center; margin: 5px;">  <p>Stemmed</p> </div> <div style="text-align: center; margin: 5px;">  <p>Inset Only</p> </div> </div>		

Screw/Pin Information

The chart below contains relevant information on various 3.2mm Pins/Screws that are compatible with the *Persona* system. If these screws/pins are used during the procedure for instrument fixation, they should be removed prior to closure as they are NOT implantable.

Screw/Pin	Screw/Pin Item #	Compatible Driver	Shipped Sterile/ Non-sterile	Quantity per Package	Single use?
	2.5mm x 2.5mm Female Hex Screw 42-5099-025-25*	 2.5mm Male Hex Driver 42-5099-025-00	Sterile	2	Yes
	7.5mm x 3.2mm Trocar Tipped Drill Pin (2.5mm hex) 00-5901-020-00		Sterile	4	Yes
	Hex Headed Screw 33mm long 00-5901-035-33	 Pin/Screw Insertor 00-5901-021-00	Sterile	2	Yes
	MIS Quad-Sparing Total Knee Headed Screw 48mm long 00-5983-040-48	 Screw Insertor/Extractor 00-5983-049-00	Sterile	1	Yes
	25mm Shorthead Holding Pin 00-5977-056-03	 Multi Pin Puller 00-5901-022-00	Non-Sterile	1	No

*Note: The 2.5mm Female Hex Screws and 2.5mm Male Hex Driver should not be used in cortical bone, as this may increase the incidence of stripping of the driver.

Optional Instrument

The Shift Block can be used to internally or externally rotate the 4-in-1 Cut Guide 2 degrees and/or shift 1mm in the anterior or posterior direction. Refer to Appendix C: Optional Instruments: Shift Block, for use.

- After final placement of the desired Anterior Referencing 4-in-1 Cut Guide, insert 3.2mm trocar-tipped pins or 3.2mm Headed Screws (see Screw Information section for examples) through the oblique holes in the Anterior Referencing 4-in-1 Cut Guide (Fig. 46).

Fig. 46

TECHNIQUE TIP

4.C

It is not recommended that the following Headed Screws are used through the oblique holes of the Anterior Referencing 4-in-1 Cut Guides, as the head of the screw may interfere with the saw blade: 00-5791-041-00, 00-5791-043-00, 00-5791-044-00, 00-5061-063-00.

Use a 1.27mm (.050 in.) thick oscillating saw blade to complete the anterior, posterior, posterior chamfer and anterior chamfer resections through the cut slots (Fig. 47). Upon completion of the cuts, use the Multi Pin Puller to remove the oblique pins.


Fig. 47

Use the Persona Slaphammer to remove the cut guide from the femur insert slaphammer rotate 1/4 turn clockwise to engage the locking feature extract (Fig. 48a & 48b).


Fig. 48a

Fig. 48b

TECHNIQUE TIP

4.D

Completing the femoral resections in the order of anterior, posterior, posterior chamfer, and then anterior chamfer, the 4-in-1 cut guide will have the greatest stability.

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