20 ottobre 2014

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**SECTION 1** | **Resect Distal Femur** | **PERSONA**  
THE PERSONALIZED KNEE SYSTEM

## Resect Distal Femur

### Assemble Adjustable Distal Resection Instrumentation

- Pull the lever on the Modular Handle and insert the IM Rod (Fig. 1). Orientation of the IM Rod will align with the polished line on the sides of the Modular Handle.

**TECHNIQUE TIP**

**1A**

*Alternately, the IM Rod may be inserted into the top of the Modular Handle to accommodate surgical preference.*

- Set the valgus angle on the Adjustable Valgus Guide by pressing the button and rotating the dial to the appropriate left or right valgus angle from 0° to 9° (Fig. 2).
- Avoid turning the Locking Knob excessively in the counterclockwise, or "unlocking", direction to prevent it from binding.
- Insert the IM Rod into the Adjustable Valgus Guide.

**Note:** If desired, align the depth markings on the IM Rod with the flat plate of the Adjustable Valgus Guide to set the IM Rod at a specific length. Inserting the IM Rod beyond the double line marking, indicated with caution symbols, may prevent assembly of the Modular Handle to the IM Rod (Fig. 3).

- Refer to Appendix C for Optional Instruments to resect the distal femur.

- Set the resection depth on the Adjustable Resection Tower by rotating the dial (Fig. 4). The '0' setting indicates a 10mm resection. Adjustments can be made in 1mm increments from 10mm to 14mm.

Fig. 1


Fig. 2

Fig. 3

Fig. 4

**Instruments**


**Note:** The '0' setting can be set to indicate a 1mm through 9mm resection depth, in 1mm increments, by assembling the corresponding Resection Plate to the Valgus Guide. **Figure 5** illustrates this with the 9mm Resection Plate. The 9mm and 8mm Resection Plates may facilitate correction of a flexion contracture and the 1mm-3mm Resection Plates can be used for recutting the distal femur, if needed.


**Fig. 5**

- Insert the Adjustable Resection Tower into the Adjustable Valgus Guide.
- Rotate the lock lever on the Adjustable Resection Tower to the unlocked, or "in-line", position and fully insert into the Cut Guide (**Fig. 6a**).

**Fig. 6a**

Flip the lock lever to the locked, or "vertical", position to secure it to the Cut Guide (**Fig. 6b**).

**Fig. 6b**

**TECHNIQUE TIP 1.B**

*The Adjustable Resection Tower is compatible with the Fixed Valgus Guide and the Fixed Resection Tower is compatible with the Adjustable Valgus Guide. These instruments can be interchanged to accommodate surgical preference.*

**Instruments**

**Establish Femoral Alignment**

- Drill the IM canal using the 8mm IM Step Drill (Fig. 7). Suction the canal to remove medullary contents.
- Insert the IM Rod and assembled distal resection instrumentation into the IM canal far enough to ensure the most accurate replication of the anatomic axis.

Fig. 7

- Set the orientation of the Adjustable Valgus Guide by placing it against the most prominent distal condyle and rotating it about the IM Rod so that the engraved lines are aligned with the epicondylar axis (Fig. 8a).
- Turn the lock knob on the Adjustable Valgus Guide clockwise, to the locked position, to secure orientation of the assembly (Fig. 8b).

**Note:** Setting rotation of the Adjustable Valgus Guide is important for creating a distal resection that matches the desired valgus angle selected. It does not set the rotation of the femoral component.


Fig. 8a

Fig. 8b

- For additional fixation, or in lieu of using the lock knob, impact the captured pin on the medial or lateral side of the Adjustable Valgus Guide until the head of the captured pin is flush with the plate (Fig. 9).


Fig. 9

**Instruments**

8mm IM Step Drill  
00-5978-014-00

Persona 8mm IM Rod  
42-5099-002-00

Persona Modular Handle  
42-5099-014-00

Persona Valgus Alignment Guide  
42-5099-004-00

Persona Adjustable Resection Tower  
42-5099-008-00

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