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
- Klein et al., 2009. 'Task Group 142 report: Quality assurance of medical accelerators.' Medical Physics, 36(nine), pp. 4197-4212.
- Bissonnette et al., 2008. 'Quality assurance for the geometric accuracy of cone-beam CT guidance in radiation therapy.' Int. J. Radiat. Oncol. Biol. Phys., 71(1), pp.S57-S61.
- Bissonnette et al., 2008. 'A quality assurance program for image quality of cone-beam CT guidance in radiation therapy.' Med. Phys., 35(5), pp.1807-1815.
- Sykes et al., 2008. 'Measurement of cone beam CT coincidence with megavoltage isocentre and image sharpness using the QUASAR™ Penta-Guide phantom.' Phys. Med. Biol., 53, pp.5275-5293.
- Cho et al., 2005. 'Accurate technique for complete geometric calibration of cone-beam computed tomography systems.' Med. Phys., 32(4), pp.968-983.

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FCO Action Notification Report

Please complete the details below and sign the appropriate acknowledgement section:

- Existing installations; Acknowledgement by the customer
- New installations: New installation confirmation by the installing Elekta or Representative employee

Please return this report to your local Elekta Office or Representative, as soon as possible and within 30 days at the latest.

***The information in this FCO has been provided to address a safety issue and therefore the customer is expected to acknowledge and accept the recommendations given, and ensure they are implemented. By refusing to implement the recommendations, the customer assumes full responsibility and liability for all matters (including costs, losses, claims, and expenses) resulting, whether directly or indirectly from not implementing such recommendations. Further the customer will hold Elekta harmless from all matters (including costs, losses, claims and expenses) resulting, whether directly or indirectly from not implementing such recommendations. Failure to sign and return the acknowledgement may affect any follow-up actions necessary for us to take, and may require Elekta to report to the Regulatory Authorities in your country.**

Classification:	Important Field Safety Notice	FCO Ref:	200 01 502 043
FCO description:	Recommended quality assurance after a calibration procedure		
Scope:	Integrity™ R1.1, R3.0, and R3.1		

Hospital:	
Device Serial No: (e.g. linac - if applicable)	Location or Site No:

Acknowledgement by customer*: This notification to be signed by the customer I acknowledge that I have read and understood this Notice and accept implementation of any given recommendations:	
Name:	Title:
Signature:	Date:

New installation confirmation: This notification to be signed by the installing Elekta or Representative employee	
Name:	Title:
Signature:	Date:

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