

Philips Healthcare

BU iXR

XCR603-140329

V.05/07

Report Form Field Safety Corrective Action

Medical Devices Vigilance System (MEDDEV 2.12/1 rev 5)

-1/4-

1. Administrative information

Destination
Name of national competent authority (NCA)

See Unit affected list

Address of national competent authority

Date of this report

Reference number assigned by the manufacturer FC072200269 / CAPA#: TW3509618

Incidence reference number and name of the co-ordinating national competent authority (if applicable)

Identify to what other national competent authorities this report was also sent See section 7

2 Information on submitter of the report

Status of submitter

- Manufacturer
- Authorised representative within EEA
- Others (identify the role):

3 Manufacturer information	
Manufacturer name Philips Medical Systems BV	
Manufacturer's contact person T. op het Veld	
Address Veenpluis 4-6, P.O. Box 10.000	
Postal code	City
5680 DA	Best
Phone	Fax
+31 6 10420792	
E-mail	Country
Twan.op.het.veld@Philips.com	The Netherlands

4 Authorised representative information				
Name of the authorised representative				
The authorised representative's contact person				
Address				
Postal code	City			
Phone	Fax			



© KONINKLIJKE PHILIPS ELECTRONICS N.V. 2008

All rights are reserved. Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical or otherwise, is prohibited without the prior written consent of the copyright owner.

Philips Healthcare

E-mail	Country					
5 National contact point information						
National contact point name						
Name of the contact person						
Address						
Postal code	City					
Phone	Fax					
E-mail	Country					
6 Medical device information						
Class AIMD Active implants						
MDD Class III	□ IVD Annex II List A					
MDD Class Ib	□ IVD Annex II List B					
MDD Class Ila	IVD Devices for self-testing					
MDD Class I	IVD General					
Nomenclature system (preferable GMDN) GMDN						
Nomenclature code 37621						
Nomenclature text						
X-Ray Imaging System Commercial name/brand name/make						
Allura Xper Systems						
Model number 722003, 722010, 722011, 722012, 722013, 722017, 722121, 722122, 722123, 722124, 722133 and 722134						
Serial number(s) and/or lot/batch number(s) See UAL						
Software version number (if applicable) N/A						
Manufacturing date/expiry date (if applicable)						
Accessories/associated device (if applicable)						
Notified body (NB) ID- number 344						

7 Description of FSCA

Background information and reason for the FSCA

We have received a customer feedback where the Monitor Ceiling Suspension system fell to its lowest position. In this specific occasion it collided with the table top. The investigation initiated concluded that the cause of the failure was an assembly error of the Actuator of the MCS.

 Description and justification of the action (corrective/preventive) The investigation initiated concluded that the cause of the failure was an assembly error of the Actuator of the MCS. The MCS must be moved –manuel or motorically- before there is a chance that it drops to its lowest position and User, patient or bystander need to stand in a close position in relation to the MCS before a hazardous situation might occur. Philips will replace the actuator of all affected systems free of charge (Mandatory FCO72200269). It is expected that this corrective action will be available as of February 2015. In the meantime, as to mitigate the risk, Philips will secure the MCS actuator of the affected systems with straps (containment action Mandatory FCO72200274). Philips representatives will contact the involved customers as soon as possible in order to schedule the implementation of this action. 						
 Advice on actions to be taken by the distributor and the user. Inform all possible system users. In order to avoid any risk for patients, users or bystanders we recommend the following until the containment action FCO72200274 or corrective action FCO72200269 has been implemented: When system in use, we recommend that the user does not: Position or move the Monitor Ceiling Suspension above the patient. Allow staff to stand under or close to the Monitor Ceiling Suspension. In addition, please avoid any unnecessary movement of the Monitor Ceiling Suspension and inform all possible System users. 						
Attached please find						
Field Safety Notion						
FSN in national la						
	pecify): Unit Affected L	ist. (UAL)				
Time schedule for the implementation of the different actions						
	in the EEA and Switze	rland are affected by t	nis FSCA			
 EEA and Switze 		_	_			
🛛 Austria	🗌 Belgium	🗌 Bulgaria	Cyprus	Czech Republic		
🛛 Denmark	Estonia	Finland	🛛 France	🖾 Germany		
Greece	🗌 Hungary	Iceland	Ireland	🛛 Italy		
Latvia	Liechtenstein	🗌 Lithuania	Luxemburg	Malta		
Netherlands	Norway	Poland	Portugal	🗌 Romania		
 □ Slovakia	☐ Slovenia	⊠ Spain	Sweden	Switzerland		
United Kingdom						
- Candidate Cour	ntries					
Croatia	Turkey					
All EEA, Candidate Countries and Switzerland.						
- Others:						
These countries outside the EEA and Switzerland are affected by this FSCA						
AU, BR, CA, CL, CN, IN, JP, MX, US, ZA.						
8 Comments						

I affirm that the information given above is correct to the best of my knowledge.

Signature:

T. op het Veld Best Name City

Date

Submission of this report does not, in itself, represent a conclusion by the manufacturer and/or authorized representative or the national competent authority that the content of this report is complete or accurate, that the medical device(s) listed failed in any manner and/or that the medical device(s) caused or contributed to the alleged death or deterioration in the state of the health of any person.