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
The screenshot shows the SPA PLUS software interface. The 'DILUTION' section is highlighted. The 'Pre Dilution Rate' field is set to 10. The 'Auto Rerun Dilution Rate High' and 'Auto Rerun Dilution Rate Low' fields are empty. The 'Diluent' is set to 99:DI1. The 'Reaction Check' is set to OFF. The 'Absorbance Limit' is set to LOW 0 and HIGH 3. The 'Monitor' section shows 0 LEVEL POINT 1 and SPAN 3. The 'Page' is 2. The 'Print', 'Prev Page', 'Next Page', 'Save', and 'Return' buttons are visible at the bottom.

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
The screenshot shows the SPA PLUS software interface. At the top, the window title is 'SPA PLUS' and the status is 'Off Line'. The 'Item Name' is '29' and 'CRP'. The 'Auto Rerun SW' section has 'On' and 'Off' radio buttons, with 'Off' selected. The 'Auto Rerun Range (Result)' section has 'Or', 'Off', 'On', and 'Off' radio buttons. Below this is a table with columns 'Lower' and 'Higher' for various sample types: Serum, Urine, Plasma, CSF, Dialysis, and Other. The 'Lower' column has a value of 0.2 for all, and the 'Higher' column has a value of 160. The 'Auto Rerun Condition (Absorbance)' section has 'Absorbance Range' with 'Lower' and 'Higher' buttons, each with 'On' and 'Off' radio buttons. The 'Prozone Range' section has 'On' and 'Off' radio buttons. The 'Bottle Size (ml)' section has two columns: '24 Items' and '36 Items', each with 'REAGENT1', 'REAGENT2 R1', and 'REAGENT2 R2' fields. At the bottom, there are buttons for 'Page : 3', 'Print', 'Prev Page', 'Next Page', 'Save', and 'Return'.

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