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Field Action Manager

Anne M. Brisson

Sr. QA Manager, Product Safety and
Performance

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Allegato 1. Confronto IPU errate (SE_018720_AE) e IPU corrette (SE_018720_AF)

<p>Incorrect</p> <p>Precautions</p> <ul style="list-style-type: none">- Do not open the package until use.- Examine the packaging for damage before using the sterile bone substitute, as they might impair sterility. This applies to both the inner primary and the outer secondary peel-off package. In removing the implant from its package, strictly observe the instructions concerning aseptic procedures.- Do not attempt to re-sterilise the unused contents of an opened package, but dispose of it. Re-sterilization alters the chronOS inject powder blend and liquid components.- Do not use chronOS inject after expiration of the use-by date printed on the package.- When preparing the chronOS for implantation, use exclusively components from one and the same package. Always mix the entire amount of powder mixture and liquid- chronOS inject must not be mixed with any additives that are not included in the package.- Always apply chronOS inject directly from the cartridge or use the delivery needle intended for the application. Never use delivery needle with a diameter inferior to 12GA. <p>Side Effects</p> <p>No Side effects due to material have been reported to date.</p> <p>Interactions</p> <p>No negative interactions have been reported to date.</p>	<p>Correct</p> <p>Precautions</p> <ul style="list-style-type: none">- An application of chronOS Inject in an enclosed defect with access to blood vessels can provoke an embolism and must be avoided.- Particular care is needed when applying chronOS Inject close to an open articular cavity. Avoid extravasations into the articular space.- Do not open the package until use.- Examine the package for damages before using the sterile bone substitute, as they might impair sterility. This applies to both the inner primary and the outer secondary peel-off package. In removing the implant from its package, strictly observe the instructions concerning aseptic procedures.- Do not attempt to re-sterilise the unused contents of an opened package, but dispose of it. Re-sterilization of the chronOS Inject Bone Void Filler powder blend and/or chronOS Inject liquid component can result in product not being sterile, and/or not meeting performance specifications and/or alters material properties.- Do not use chronOS Inject Bone Void Filler after expiration of the use-by date printed on the package.- When preparing chronOS Inject Bone Void Filler for implantation, use exclusively components from one and the same package. Always mix the entire amount of powder mixture and liquid.- chronOS Inject Bone Void Filler must not be mixed with any additives that are not included in the package.- Always apply chronOS Inject Bone Void Filler directly from the cartridge and use the delivery needle intended for the application. Never use delivery needles with a diameter inferior to 12 ga (gauge).	<p>Correct</p> <p>Possible Side Effects, Possible Adverse Effects and Potential Complications</p> <ul style="list-style-type: none">- Non union or delayed union, which may lead to a failure of the implant- Pain, discomfort, abnormal sensation, or palpability due to the presence of the device- Increased fibrous tissue response around fracture site and/or the implant <p>Apart from these possible adverse effects there is also the risk of complications associated with any surgical procedure featuring bony defects such as, but not limited to, necrosis of bone, infection, nerve damage and pain which may not be related to the implant.</p> <p>In general, good tissue response of tricalcium phosphate/brushite (dicalcium phosphate dihydrate) implants in bone is supported by experimental and clinical data.</p> <p>Nevertheless, the following complications are possible:</p> <ul style="list-style-type: none">- Fragment displacement as a result of use in inappropriate indications- Neurovascular injuries caused by surgical trauma- Foreign body reactions- Allergic reactions- Inflammatory reactions- Infections that can lead to failure of the procedure- General complications caused by invasive surgery
---	---	--

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Synthes GmbH