v.01.13

1. Administrative Information		
Destination		
Ministry of Health		
Vigilance on Medical Devices Via Giorgio Ribotta 5, IT - 00144 Roma		
Italy		
Type of Report		
☐ Initial report☐ Follow up report		
Final report		
Date of this Report		
30 October 2015		
Reference Number Assigned by Manufacture	,	
21833502-01/19/2015-001-R		
FSCA Reference Number Assigned by NCA		
2015/001/029/071/001		
Incidence Reference Number Assigned by NC	^^	
N/A	<b>/A</b>	
Name of the Coordinating National Competer	t Authority (if applicable)	
MHRA	, ,	
2. Information on Submitter of the Report		
Status of submitter		
☐ Manufacturer		
Authorised Representative within EEA, \$	Switzerland and Turkey	
Others: (identify the role):		
Smiths Medical Risk Management Specialist on I	Manufacturer's Behalf	
3. Manufacturer Information		
Manufacturer Name		
Smiths Medical ASD, Inc.		
Manufacturer's Contact Person		
Tim Giguere		
Address		
1265 Grey Fox Road  Postal Code	City	
55112	City St. Paul	
Phone	Fax	
651 628 7477	n/a	
E-mail	Country	
tim.giguere@smiths-medical.com	USA	
4. Authorised Representative Information		
Name of Authorised Representative		
Smiths Medical International Ltd.		

v.01.13

Authorised Representative's Contact Person			
Marco Savino			
Address 1500 Eureka Park, Lower Pemberton, Kent			
Postal Code	Postal Code		
TN25 4BF	TN25 4BF		
Phone	Phone		
39 0773 4084810	39 0773 4084804		
E-mail	Country		
eu.rep@smiths-medical.com	United Kingdom		
5. National Contact Point Information			
National Contact Point Name			
Smiths Medical International Ltd.			
Name of the Contact Person			
Marco Savino			
Address			
1500 Eureka Park, Lower Pemberton, Kent			
Postal Code	Postal Code		
TN25 4BF	TN25 4BF		
Phone	Phone		
39 0773 4084810	39 0773 4084804		
E-mail	Country		
eu.rep@smiths-medical.com United Kingdom			
6. Medical Device Information			
Class			
☐ AIMD Active Implants			
☐ MDD Class III ☐ IVD Annex II Li:	st A		
☐ MDD Class IIb ☐ IVD Annex II Lis	st B		
	Self-Testing		
☐ MDD Class I ☐ IVD General	-		
Nomenclature System (preferable GMDN)  Nomenclature Code			
GMDN	35127		
Nomenclature Text			
Intravenous fluid container, single use			
Commercial Name/ Brand Name / Make			
CADD™ Medication Cassette Reservoir			
Model Number	Catalogue Number		
N/A	21-7001-24, 21-7301-24		
Serial Number(s)	Lot/ Batch Number(s)		
N/A	21-7001-24, Lot Numbers 14X-297 and 14X-323.		
Davisa Manufacturing Data	21-7301-24, Lot Number 14X-324		
<b>Device Manufacturing Date</b> 14X297 12-June-2014	Expiry Date June 2019		
14A237 12-JUNE-2014	Julie 2019		

v.01.13

14X323 21-June-2014	
14X324 21-June-2014	

## **Software Version Number (if applicable)**

N/A

#### Accessories/ Associated Device (if applicable)

ΝΙ/Δ

#### Notified Body (NB) ID Number

0473

#### 7. Description of FSCA

#### **Background Information and Reason for the FSCA:**

Smiths Medical has become aware of an issue with specific lots of 50mL CADD™ Medication Cassette Reservoirs ("Cassette"). Some Cassettes may leak at the sealing area of the pump tube and medication bag. Smiths Medical has received no reports of serious injury or death related to this issue.

Examination of Cassettes returned for investigation confirmed leakage at the sealing area of the pump tube and medication bag. Investigation found that the complaints with confirmed leaking were limited to 3 finished goods lots. Despite an in-depth investigation, we have not identified a definitive root cause. Bag and tube material, set-up related training and technician training were all investigated and eliminated as potential root causes. The most probable root causes found during the investigation are machine and process related:

- 1. Damaged plate holder mandrels;
- 2. Failures in the grounding mechanism of the power supplied to the dies; and/ or
- 3. Worn motor brushes.

#### **Description and Justification of the Action (Corrective/ Preventive):**

As discussed within the Risk Analysis Summary submitted with the original notification; despite indepth investigation, the definite root cause for the leakage was not identified. However, the root cause was isolated to manufacturing. The issue occurred during production when the pump tube and medication bag were sealed together (this "sealing area" of the product is circled in photo below).



The potential root causes were defined as:

- 1. Damaged plate holder mandrels (equipment related);
- 2. Failures in the grounding mechanism of the power supplied to the dies (equipment and process related).

In order to address these 2 potential root causes identified, Smiths Medical has implemented the following actions:

Preventive Maintenance Procedure for the bag machine was updated to include weekly inspection of the plate holder mandrels and ground strap. The ground strap is the mechanism used to supply power to the dies during the sealing operation.

v.01.13

The actions taken address potential root causes 1 and 2 above.			
Advice on Actions to be Taken by the Distributor and the User:  All consignees were sent an Urgent Field Safety Notice via mail service to notify them of this Field Action. Consignees were instructed to return the product for credit or replacement. Distributors were instructed to notify their customers. The Urgent Field Safety Notice included a Confirmation Form that consignees were instructed to send back to Smiths Medical for carrying out the action and tracking effectiveness.			
Progress of FSCA with Reconciliation Data (Mar	ndatory for a Final FSCA)		
All affected consignees in Italy have sent Smith's Medical their acknowledgement of receipt of this notice and returned affected stock. This action has been completed for Italy.			
Attached Please Find  Field Safety Notice (FSN) in English FSN in National Language Others (please specify)	SN Status Draft Final		
Time Schedule for the Implementation of the Dif	ferent Actions:		
Please see above.			
These Countries Within the EEA, Switzerland, and Turkey are Affected by this FSCA:			
Within the EEA, Switzerland, and Turkey:  △AT △BE □BG △CH □CY △CZ △DE □  △FI △FR △GB □GR △HU △IE □IS □  □LU □LV □MT △NL □NO □PL △PT  □SK □TR	⊠IT □LI □LT		
Candidate Countries:			
□HR			
☐ ALL EEA, Candidate Countries, Switzerland, and Turkey			
Others:			
AE, AU, CA, CO, ID, SG, US, ZA			
8. Comments: As Smiths Medical's Notified Body, Intertek, is located in the UK, Smiths Medical recognizes the MHRA as the Lead Competent Authority for this product.			
I affirm that the information given above is correct to the best of my knowledge.			
Ale Riebe	00.0 1 1 0045		
Signature	30 October 2015 Date		
Ellen Riebe	St. Paul, MN		
Name	City		

v.01.13

Submission of this report does not, in itself, represent a conclusion by the manufacturer and / or authorized representative or the National Competent Authority that the content of this report is complete or accurate, that the medical device(s) listed failed in any manner and/or that the medical device(s) caused or contributed to the alleged death or deterioration in the state of the health of any person.