
CORRECTIVE AND PREVENTIVE ACTION FORM

No: HT20160425

Customer: Orange Medical
PI number: HTF150212IT
Shipping date: May 16th 2015
Date receiving complaint from customer: 18th April 2016

| | |
|---|---|
| Products | Endotracheal tube cuffed |
| size | 8.0 mm |
| Complaints from customer | one hospital has reported following malfunctioning on tracheal tube # 111-080: "After the laryngoscopy the device was positioned inside the trachea, the balloon was inflated with the pressure gauge, but there were several leakages in the ventilation circuit. The device was replaced with another device (same type). When the tracheal tube was pulled out, the "warning balloon" remained inflated but the end balloon that should keep the tube positioned, was completely deflated" |
| Root Cause | After checking sample sent from Ms Simona Schirato Orange Medical, we found that the "Inflation line" was completely squashed . Please see picture attached. This caused the "warning balloon" can be inflated and the "end balloon" was deflated. Issued by: Huang yuanhua, QC date: 2016-04-25 |
| Cause Analysis | 1. If the inflation line was founded squashed when opening the banana pouch, then, the reason should be "the inflation line was squashed during automatic packing process. And our people did not find the defected product during "light inspection". 2. If the inflation line was squashed in the operation room, the "warning balloon" can be inflated, but the "end balloon" can not be inflated or deflated. |
| Corrections | Corrections for Cause analysis 1: a. To be more careful during position ET tube into banana pouch for automatic packing b. Asking Light Inspection staff to add another checking on inflation line, before packing into middle boxes, so as to make sure 100% of inflation lines are not squashed Issued by: Liao quangen, QC manager Date: 2016-04-25 |
| Suggestions | 1. Usually, before intubation, doctors will test the balloon by inflating and deflating. If the inflation line was squashed or cannot work, they will change for another tube. 2. Please kindly do us a favor for checking when the inflation line was squashed, before operation or during. |
| Approved: Sheng Wenping  date: 2016-04-25 | |

Picture of ET tube sample sent from Ms Simona Schirato Orange Medical

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