LabStyle Innovation Ltd.

Email Subject: Data Presented in Dario's Android Application LOGBOOK May Potentially be Incorrect

Important field safety notice

Data Presented in Dario's Android Application LOGBOOK May Potentially be Incorrect

Dear User,

We would like to inform you about a potential issue we have identified with Android App versions of Dario - Blood Glucose Tracker (4.3.0-4.3.2), regarding duplicate logging of a blood glucose level reading from Dario™ Blood Glucose Monitoring System for Android.

Please note that this potential issue only affects users of Dario - Blood Glucose Tracker Android App versions **4.3.0-4.3.2** who are using the Dario™ Blood Glucose Monitoring System for Android.

You are not affected and can use Dario's system without any limitations, in case you are using:

- An iPhone with the Dario™ Blood Glucose Monitoring System or the Dario™ LC Blood Glucose Monitoring System
- Or using any other Android App versions.

Please note that the blood glucose level reading seen in the measurement screen (inside the circle), while taking a measurement, is unaffected and displayed correctly.

The issue may cause blood glucose level readings that are transferred from the measurement screen to the 'LOGBOOK' may be replaced by an older blood glucose measurement and be displayed as a duplicate reading.

If inaccurate blood glucose level reading from the 'DATA ENTRY' or 'LOGBOOK' screens are used to determine the insulin dose or make treatment decisions, this may lead to incorrect dosing decisions. As safety is of utmost importance to us, LabStyle Innovations Ltd. is informing users of Dario™ Blood Glucose Monitoring System for Android about this potential issue and is updating the Dario - Blood Glucose Tracker App to resolve this problem.

What do you need to do?

Please update the latest Dario - Blood Glucose Tracker Android App version as follows -

- Open the Google Play store
- Search for Dario Blood Glucose Tracker App
- Press the 'Uninstall' button
- Press 'OK' to confirm the question "Do you want to uninstall this app?"
- Once completed, press 'Install' and verify you installed **4.3.3** version (or above) by pressing the 'About this App' link.

We would like to thank you for your co-operation and apologize for any inconvenience this may cause. If you have any questions, please feel free to contact our customer support. Details can be found at: <http://mydario.com>.

Best Regards,

Michal HersHKovitz

Michal HersHKovitz

VP Clinical, Regulatory Affairs & Quality Assurance

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