FSN Ref: M18022

FSCA Ref: M18022

Date: 23 October 2019

**Urgent Field Safety Notice**  
**Clinodigit EVO**

**For Attention of:** Dealers and Users with Cinodigit EVO


Contact details of local representative (name, e-mail, telephone, address etc.)*
This could be a distributor or local branch of the manufacturer. To be added at the appropriate stage in the different local languages

**Urgent Field Safety Notice (FSN)**  
**Clinodigit EVO**  
**Angulation chain issue**

<b>1. Information on Affected Devices*</b>	
1.	1. <b>Device Type(s)*</b> Remote controlled table
1.	2. <b>Commercial name(s)</b> Clinodigit EVO
1.	3. <b>Unique Device Identifier(s) (UDI-DI)</b> N/A
1.	4. <b>Primary clinical purpose of device(s)*</b> Medical purpose is the diagnosis, visualization and documentation of various human diseases or injuries. The system permit the following Clinical applications: Routine X-ray examinations and special examinations on patients in the standing, seated or lying positions. This includes standard radiography and standard fluoroscopy procedures: Lung fluoroscopy, Barium studies of the upper and lower gastro intestinal tract, Other examinations of the digestive tract, Urinary studies (IVP, Cystourethrography), Myelography, Arthrography, Hysterosalpingography, Phlebography, Examinations of the gall-bladder / ERCP, Pediatric examinations, Tomography, X-ray examinations of the skeleton (including skull, ribs, spine, pelvis, upper extremities, lower extremities, etc.), X-ray examinations of the lung, X-ray examinations of soft tissue such as abdomen, Stitching examinations (optional - mostly for lower limbs and spine).
1.	5. <b>Device Model/Catalogue/part number(s)*</b> 9784520013
1.	6. <b>Software version</b> not relevant
1.	7. <b>Affected serial or lot number range</b> 15101490, 15111497, 16011509
1.	8. <b>Associated devices</b> N/A

<b>2 Reason for Field Safety Corrective Action (FSCA)*</b>	
2.	1. <b>Description of the product problem*</b> It was found that the angulation chain is very close to the design limit for metal fatigue and, under severe workload conditions might fail before the expected lifetime of 10 years. Therefore, the Ifu and the Service Manual need to be amended to prescribe the replacement of the angulation chain after 8 years of service life of the unit.
2.	2. <b>Hazard giving rise to the FSCA*</b> In case of failure of the chain, the column falls and stops hitting the shoulders of the DRF carrier / SFD as shown in the following picture, so the column does not come in contact with the tabletop since the distance between the tabletop and the collimator, in this condition, is about 0.5 m:

	 <p>Generally, the patient lays on the tabletop and the operator is away from the table, or in the front side.</p> <p>In case this advice is followed, the residual risk is lowered to the initially estimated level, with likelihood lowered to “Improbable” level.</p>
2.	<p><b>3. Probability of problem arising</b></p>
	<p>The likelihood is considered Remote (only under severe workload conditions)</p>
2.	<p><b>4. Predicted risk to patient/users</b></p>
	<p>The identified hazard is: Mechanical hazard (collision, squeezing) due to the column fall due to a broken chain due to fatigue. The risk has severity “serious” and probability “remote”</p>
2.	<p><b>5. Further information to help characterise the problem</b></p>
	<p>Only one case happened, without involvement of patient or operator</p>
2.	<p><b>6. Background on Issue</b></p>
	<p>Following information received from field, it was found that the number of tilting cycles for some units used in high-workload environments had to be increased with respect to the initial specifications. After laboratory testing of the current angulation chain, it was found that it is very close to the design limit for metal fatigue and, under severe workload conditions might fail before the expected lifetime of 10 years</p>
2.	<p><b>7. Other information relevant to FSCA</b></p>
	<p>N/A</p>

	<p><b>3. Type of Action to mitigate the risk*</b></p>	
3.	<p><b>1. Action To Be Taken by the User*</b></p>	
	<p>Take note of amendment/reinforcement of Instructions For Use (IFU) and Service Manual, attached to this FSN. Replace the attached pages in the User Manual and in the Service Manual, prescribing the periodic replacement of the angulation chain. Fill the table present in the addendum; In the cell “Next replacement date” write the date adding 8 years to the replacement date</p>	
3.	<p><b>2. By when should the action be completed?</b></p>	<p>as soon as possible</p>

3.	<b>3.</b> Particular considerations for Diagnostic Imaging device Is follow-up of patients or review of patients' previous results recommended? No The issue has no impact on the images taken	
3.	<b>4.</b> Is customer Reply Required? * (If yes, form attached specifying deadline for return)	Yes
3.	<b>5. Action Being Taken by the Manufacturer</b> An ECO has been implemented in order to adopt a new type of angulation chain, whose design limit for metal fatigue exceeds the expected lifetime of 10 years. Instructions For Use (IFU) and Service Manual have been updated to prescribe the replacement of the angulation chain after 10 years of service life of the unit.	
3	<b>6.</b> By when should the action be completed?	The action is already completed: first unit with new type of chain and updated IfU and Service Manual is the unit with s/n 19102065 manufactured on October 2019
3.	<b>7.</b> Is the FSN required to be communicated to the patient /lay user?	No

<b>4. General Information*</b>		
4.	<b>1. FSN Type*</b>	New
4.	2. For updated FSN, reference number and date of previous FSN	N.A.
4.	3. For Updated FSN, key new information as follows:	N.A.
4.	<b>4. Further advice or information already expected in follow-up FSN? *</b>	No
4.	5. If follow-up FSN expected, what is the further advice expected to relate to:	N.A.
4.	6. Anticipated timescale for follow-up FSN	N.A.
4.	<b>7. Manufacturer information</b> (For contact details of local representative refer to page 1 of this FSN)	
	a. Company Name	<b>Villa Sistemi Medicali S.p.A.</b>
	b. Address	<b>via delle Azalee 3 - 20090 Buccinasco (MI) - ITALY</b>
4.	<b>8. The Competent (Regulatory) Authority of your country has been informed about this communication to customers. *</b>	
4.	9. List of attachments/appendices:	1. IfU and Service Manual addendum 2. Distributor and User reply form
4.	10. Name/Signature	<b>Paolo Casagrande Santin</b> <b>Quality Assurance Manager</b> 

FSN Ref: M18022

FSCA Ref: M18022

	<b>Transmission of this Field Safety Notice</b>
	<p>This notice needs to be passed on all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred.</p> <p>Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.</p> <p>Please report all device-related incidents to the manufacturer, distributor or local representative, and the national Competent Authority if appropriate, as this provides important feedback.</p>

FSN Ref: M18022

FSCA Ref: M18022

## **Attachment 1**

### **IfU and Service Manual Addendum**

FSN Ref: M18022

FSCA Ref: M18022

## **Addendum to**

# **Clinodigit EVO**

## **User's Manual**

The User's Manual provided with your equipment does not include all the Maintenance information for the device.

Please insert the following pages at the beginning of chapter Maintenance of your User's Manual.

All other specifications are unchanged with respect to the ones provided in your Manual.

## User Manual Addendum

The expected service life for Clinodigit EVO angulation chain is 8 years.  
 Contact your service organization in order to plan the intervention once the unit is close to reach the next replacement date given in the table below.  
 At replacement, FSE will fill the column “Next replacement date” adding 8 years to the replacement date in both User and Service Manual

Angulation chain replacement date (dd – mm – yyyy)	Reason of the replacement  corrective maintenance / planned maintenance	next replacement date mm/yyyy	Executed by Organization (FSE Name)	Date / Signature

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# **Clinodigit EVO**

## **Service Manuals**

The Service Manual provided with your equipment does not include all the Maintenance information for the device.

Please insert the following pages at the beginning of chapter “Maintenance and Cleaning” of your Service Manual.

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## Service Manual Addendum

### 9.2 Parts to be replaced periodically

The expected lifetime of Clinodigit EVO angulation chain is 8 years. After the operation time of 8 years, the following part must be replaced:

Code	Description	Replace after	Installation date
6690702200	Angulation chain	Every 8 years	

At replacement, FSE have to fill the column “Next replacement date” adding 8 years to the replacement date in both User and Service Manual

Angulation chain replacement date (dd - mm - yyyy)	Reason of the replacement  corrective maintenance / planned maintenance	next replacement date mm/yyyy	Executed by Organization (FSE Name)	Date / Signature

FSN Ref: M18022

FSCA Ref: M18022

## **Attachment 2**

### **Distributor and User reply forms**

FSN Ref: M18022

FSCA Ref: M18022

## Distributor/Importer Reply Form

<b>1. Field Safety Notice (FSN) information</b>	
FSN Reference number*	M18022
FSN Date*	23 October 2019
Product/ Device name*	Clinodigit EVO
Product Code(s)	
Batch/Serial Number (s)	

<b>2. Distributor/Importer Details</b>	
Company Name*	
Account Number	
Address*	
Contact Name*	
Title or Function	
Telephone number*	
Email*	

<b>3. Return acknowledgement to Sender</b>	
Email	<a href="mailto:p.casagrande@villasm.com">p.casagrande@villasm.com</a> ; <a href="mailto:vsmservice@villasm.com">vsmservice@villasm.com</a>
Postal Address	Villa Sistemi Medicali S.p.A. via delle Azalee 3, 20090 Buccinasco (MI), ITALY  to the attention of Paolo Casagrande Santin
Deadline for returning the Distributor/Importer reply form*	31 December 2019

<b>4. Distributors/Importers (Tick all that apply)</b>		
<input type="checkbox"/>	*I confirm the receipt, the reading and understanding of the Field Safety Notice.	Distributor/Importer to complete or enter N/A
<input type="checkbox"/>	I have checked my stock and quarantined inventory	Distributor/Importer to enter quantity and date
<input type="checkbox"/>	I have identified customers that received or may have received this device	
<input type="checkbox"/>	I have attached customer list	
<input type="checkbox"/>	I have informed the identified customers of this FSN	Date of communication:
<input type="checkbox"/>	I have received confirmation of reply from all identified customers	
<input type="checkbox"/>	Neither I nor any of my customers has any affected devices in inventory	
Print Name*		Distributor/Importer print name here
Signature*		Distributor/Importer sign Here
Date *		

Mandatory fields are marked with \*

It is important that your organisation takes the actions detailed in the FSN and confirms that you have received the FSN.

Your organisation's reply is the evidence we need to monitor the progress of the corrective actions.

FSN Ref: M18022

FSCA Ref: M18022

## Customer Reply Form

1. Field Safety Notice (FSN) information	
FSN Reference number*	M18022
FSN Date*	23 October 2019
Product/ Device name*	Clinodigit EVO
Product Code(s)	
Batch/Serial Number (s)	

2. Customer Details	
Account Number	
Healthcare Organisation Name*	
Organisation Address*	
Department/Unit	
Contact Name*	
Title or Function	
Telephone number*	
Email*	

3. Customer action undertaken on behalf of Healthcare Organisation		
<input type="checkbox"/>	I confirm receipt of the Field Safety Notice and that I read and understood its content.	Customer to complete or enter N/A
<input type="checkbox"/>	I performed all actions requested by the FSN.	Customer to complete or enter N/A
<input type="checkbox"/>	The information and required actions have been brought to the attention of all relevant users and executed.	Customer to complete or enter N/A
<input type="checkbox"/>	Other Action (Define):	
<input type="checkbox"/>	I do not have any affected devices.	Customer to complete or enter N/A

FSN Ref: M18022

FSCA Ref: M18022

<input type="checkbox"/>	I have a query please contact me (e.g. need for replacement of the product).	Customer to enter contact details if different from above and brief description of query
Print Name*		Customer print name here
Signature*		Customer sign here
Date*		

<b>4. Return acknowledgement to sender</b>	
Email	<a href="mailto:p.casagrande@villasm.com">p.casagrande@villasm.com</a> ; <a href="mailto:vsmervice@villasm.com">vsmervice@villasm.com</a>
Postal Address	Villa Sistemi Medicali S.p.A. via delle Azalee 3, 20090 Buccinasco (MI), ITALY to the attention of Paolo Casagrande Santin
Fax	+39 02 48859 303; +39 02 48859 222
Deadline for returning the customer reply form*	31 December 2019

Mandatory fields are marked with \*

It is important that your organisation takes the actions detailed in the FSN and confirms that you have received the FSN.

Your organisation's reply is the evidence we need to monitor the progress of the corrective actions.