

1  
2  
3  
4  
5  
6  
7  
8  
9  
10  
11  
12  
13  
14  
15  
16  
17  
18  
19  
20  
21  
22  
23  
24



*Ministero della Salute*



*Agenzia Nazionale per i Servizi Sanitari Regionali*

**Documento in consultazione pubblica**

Horizon Scanning report No.27

**SIRIO H3: virtual navigation system to support radiological procedures with percutaneous access**

**JUNE 2020**

25 **Methods**

26 Agenas is a public body. Its mission is to promote innovation and development within the Italian national  
 27 healthcare service and provide an Early Awareness and Alert (EAA) service by Horizon Scanning (HS) activities  
 28 in the field of new and emerging health technologies. A full description of the methods used for the production  
 29 of the present HS report can be found at [www.agenas.it](http://www.agenas.it)

30 *This document should be cited as follow:*

31 Chiarolla E, Orzella L, Amicosante AMV, Liotta A, Barbero S. SIRIO H3: virtual navigation system to support  
 32 radiological procedures with percutaneous access. Agenas, Agenzia nazionale per i servizi sanitari regionali.  
 33 Rome, June 2020

34

35 *For further information contained in this report please contact:*

36 Agenas – Agenzia nazionale per i servizi sanitari regionali  
 37 Ufficio HTA: Innovazione e sviluppo a supporto delle Regioni  
 38 Via Piemonte, 60 – 00187 Roma  
 39 e-mail: [hta@agenas.it](mailto:hta@agenas.it)

40

41 **Funding**

42 The production of the present report was made possible by financial contributions from the Italian Ministry of  
 43 Health and Agenas. Agenas takes sole responsibility for the final form and content of this report. The views  
 44 expressed herein do not necessarily represent the views of the Italian Ministry of Health or any regional  
 45 government. The intellectual content of the report is property of Agenas.

46

47 **Limitations**

48 This report is based on information available when the research was done and does not contain data on  
 49 subsequent developments or improvements of the evaluated technology. The observations made on  
 50 effectiveness, safety or cost-effectiveness of the technology evaluated in the report are to be considered  
 51 current, but may change as more evidence becomes available if an update of the document is commissioned.

52

53 **Authors**

54 *This HS report was prepared by:*

55 Emilio Chiarolla (Agenas)  
 56 Letizia Orzella (Agenas)  
 57 Anna Maria Vincenza Amicosante (Agenas)  
 58 Angela Liotta (Agenas)  
 59 Stefano Barbero (AON di Alessandria)

60 **Declaration of Conflict of Interest**

61 The authors declare that they will not receive either benefits or harm from the publication of this report.  
 62 None of the authors have or have held shares, consultancies or personal relationships with any of the  
 63 producers of the devices assessed in this document.

64

65 **Acknowledgements**

66 **Internal reviewer:** Simona Paone (Agenas)

67 **External reviewer:** Bruno Beomonte Zobel (Università Campus Bio-Medico di Roma, Facoltà di Medicina e  
 68 Chirurgia)

69 **Industry:** Walter Contento (Masmec SpA)

70

71

72

73

74  
75**HORIZON SCANNING REPORT – No. 27****Name of the technology/procedure:****SIRIO H3: virtual navigation system to support radiological procedures with percutaneous access**76 **Target population**

77 Patient candidates to tissue biopsy or drug/ablation therapy in a specific pathologic area are the  
78 population that could benefit from the SIRIO H3.

79 **Description of the procedure and technology**

80 X-ray, ultrasound, Computed Tomography (CT), Magnetic Resonance Imaging (MRI) scanning are  
81 often used as a guide for mini invasive procedures through the percutaneous insertion of a biopsy  
82 needle [ACR, 2018]. Biopsy is a medical procedure consisting of taking a small sample of tissue or  
83 cells for further microscopic evaluation. Transthoracic needle biopsy (TNB) is a minimally invasive  
84 procedure useful when imaging tests cannot exclude the malignancy of a nodule lung and/or it is  
85 difficult to reach the lesion by bronchoscopy or other methods [ACR, 2018]. Needle biopsy is less  
86 invasive than surgical ones and may not require general anesthesia.

87 The standard technique is the introduction of a needle for CT-guided biopsy during CT scan to allow  
88 the exact location of the lesion. However, this technique has some limitations with regards to high  
89 adsorbed doses, execution time and the inability to visualize in real time needle progression in  
90 reaching the lesion [Iannelli, 2018]. Virtual navigation systems are emerging tools in percutaneous  
91 imaging-guided procedures using electro-magnetic, optical or hybrid tracking during operations  
92 [Faiella, 2018]. Virtual CT-guided navigation system for biopsy provides a 3D reconstruction from a  
93 data set of acquired CT images through automatic procedures.

94 In comparison with traditional techniques, the use of CT-guided navigation system technologies  
95 allows faster operations, with greater efficacy and also with a reduction in radiation dose for both to  
96 patient and medical staff.

97 SIRIO H3 is a CT-guided navigation system designed to support radiological interventional  
98 procedures. Moreover, it is an aid for percutaneous interventions such as biopsies and/or thermal  
99 ablations in different anatomical areas, such as the lungs, bones and kidneys.

100 **Clinical importance and burden of disease**

101 Biopsies are performed to diagnose or to confirm infections, inflammations and tumors affecting  
102 various internal organs, as well as to evaluate the course of pathologies and to establish the possible  
103 therapies to be administered. It is commonly used to diagnose liver diseases, acute or chronic kidney  
104 diseases, pathologies affecting the lungs and the pleura, the prostate etc. Biopsies are mostly  
105 performed for insight into possible cancerous conditions.

106 It was estimated that around 371,000 new cancer diagnoses (excluding skin cancer) occurred in  
107 Italy in 2019 (53% among men and 47% among women). Breast cancer, lung carcinomas in women  
108 and, in both genders, pancreas, thyroid and melanomas (especially in southern Italy) are growing.

109 In the entire population, with the exception of skin carcinomas, the five most frequent tumor sites  
110 involve breast (14%), colorectal (13%), lung (11%), prostate (10%) and bladder (8%) [AIRTUM,  
111 2019]. In 2017, deaths due to tumors were 179,502 (100,003 among men and 79,499 among  
112 women) Cancer conditions with the highest number of deaths were lung cancer (33,759), followed

113 by colorectal cancer (19,355), breast (12,942), pancreas (12,347) and liver (9,214 9214) [ISTAT,  
114 2017].

115 In the period 2005-2009, compared to previous one (2000-2004), the proportion of survivors  
116 increased for both men (54% vs 51%) and women (63% vs 60%) since initial cancer diagnosis. In  
117 2019, lung cancer was the second most frequent condition in men (15%) and the third most frequent  
118 in women (12%). Estimates of new lung cancer diagnosis is around to 42,500 cases (29,500 in men  
119 and 13,000 in women). Lung cancer is the leading cause of death (12%) of all malignancies in the  
120 whole Italian population and the 5-year survival rate of lung cancer patients is equal to 16%.  
121 [AIRTUM, 2019].

122 Early diagnosis of a suspected lesion allows treatment to be started previously and can improve the  
123 prognosis in terms of survival and quality of life.

## 124 **Products, manufacturers, distributors and approval**

125 SIRIO H3 by Masmec SpA received the CE mark in 2013 and the FDA 510(k) clearance in November  
126 2017. In Italy, according to the *Classificazione Nazionale dei Dispositivi Medici (CND)* SIRIO H3 is  
127 classified under the class "CND: Z12011401 – SURGICAL NAVIGATION SYSTEM" and registered with  
128 the number 358838 within the Italian National medical device database (BD/RDM). SIRIO H3 is  
129 classified in I risk class.

130 SIRIO H3 is a medical device supporting minimally invasive interventional radiology procedures  
131 (under CT guidance) performed on the thoracic district and that involves the insertion of a probe  
132 needle (SIRIO H3 IFU). In particular, the system acts as an aid for biopsy samplings and  
133 thermoablation procedures of pulmonary nodules. SIRIO H3 is indicated in those cases defined by  
134 medical staff in which clinical conditions require an interventional radiology procedure in the thoracic  
135 district through the insertion of a probe needle (SIRIO H3 IFU).

136 According to the manufacturer, SIRIO H3 is based on 3D reconstructed CT images, and allows  
137 navigation of the anatomical district of interest, to identify the needle trajectory and to track it until  
138 the neoformation/lesion is reached, even for small lesions (less than 1cm).

139 SIRIO system is composed of an instrumented column, an intervention kit, an infrared optical sensor,  
140 a disposable sterile kit, a posture tracking system, breathing sensors and a needle support [Caparelli,  
141 2015]. The instrumented column is equipped with a visualization and image processing unit, where  
142 the reconstructed 3D model of the anatomical district of interest (virtual space) is visualized and on  
143 which the procedure is performed, and an infrared camera allowing to follow the intervention needle  
144 and the patient in the navigable virtual space using a tracking technology. The infrared camera  
145 receives the emitted radiation by the reflective elements placed on the intervention kit and calculates  
146 their position in space. Starting from tomographic scan DICOM images, the system reconstructs the  
147 three-dimensional model of the anatomical district of interest. The infrared system detects the  
148 relative position of the needle and the advancement inside the patient's chest through a real-time  
149 tracking system [Caparelli, 2015]. The virtual model helps the interventional radiologist to assess  
150 the area and the initial insertion trajectory with better approximation, as well as any corrections-

151 The tracking posture procedure monitors the correspondence between reality and three-dimensional  
152 reconstruction using single-use passive markers placed in appropriate areas of the patient. The  
153 accuracy of path in that system is <2 mm [Caparelli, 2015]. The virtual navigation is intended to  
154 drastically reduce, time procedure, number of CT scans and consequently adsorbed radiation dose  
155 by patient and healthcare personnel.  
156

Product name [Manufacturer]	Distributor	CE Mark	BD/RDM	FDA
SIRIO H3 [Masmec SpA]	Masmec SpA	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

## 157 **Roll out in Italy**

158 SIRIO's first version was launched in Italy in 2010 and a newer version in 2013. According to the  
 159 manufacturer, in Italy SIRIO H3 is used in Piemonte, Lombardia, Toscana, Emilia-Romagna, Abruzzo,  
 160 Lazio, Campania, Puglia e Basilicata, mainly in public hospitals. Only two healthcare Centres are  
 161 private.  
 162

<input type="checkbox"/> Pre-marketing	<input type="checkbox"/> On the market for 1-6 months	<input type="checkbox"/> On the market for 7-12 months
<input checked="" type="checkbox"/> On the market for more than 12 months	<input type="checkbox"/> Not identified	

163 **Setting**

164 SIRIO H3 is used in a CT interventional radiology room.  
 165

<input type="checkbox"/> Home	<input checked="" type="checkbox"/> Hospital	<input type="checkbox"/> Outpatient
<input type="checkbox"/> Accident and Emergency	<input type="checkbox"/> Other:	

166 **Comparators**

167 The standard of care is represented by procedure using only CT as aid for minimally invasive surgery.  
 168

169 **Effectiveness and safety**

170 We searched the major databases including MEDLINE (search date 29/04/2020), Clinicaltrials.gov  
 171 and Cochrane Library looking for studies on humans, published from 2010 to 2020, in Italian or  
 172 English, reporting on effectiveness and safety of SIRIO H3 MASMEC for all those who have to  
 173 undergo a biopsy or need drug/ablation therapy in a defined area. Search results (n=59: 43 Pubmed,  
 174 13 Cochrane, 3 Other sources, 1 citation duplication) were screened by reading title and abstract.  
 175 Twelve articles were considered eligible for full-text analysis. No clinical trials or systematic reviews  
 176 were identified. Six studies were pertinent to SIRIO. The list of studies read in full text, included and  
 177 excluded (with reasons for exclusion) is in Appendix 1.

178 *Comparative studies*

179 The article of Grasso et al [Grasso, 2012] reports a prospective study in performing percutaneous  
 180 lung biopsy (PLB) with the enrollment of 180 consecutive subjects, over a period of 14 months.  
 181 Patients, not randomized, were equally divided into a control group (group C) and a study group  
 182 (group S). Parameters were measured in both groups of patients as follows: a) time of execution of  
 183 the procedure; b) total number of CT scans performed and x-ray dose; c) diameter of the lesion  
 184 biopsied; d) technical success (achievement of the target lesion) and diagnostic success (final  
 185 pathological diagnosis); e) procedural complications. The time of execution the procedure has shown  
 186 to be substantially shorter than that of procedures performed with the traditional technique (mean  
 187 execution time 14 minutes vs 24 minutes), with a statistically significant difference (P<0.001). The  
 188 number of CT scans performed (z = 5.64) and, consequently, the mean x-ray dose (CTDIvol)  
 189 administered to patients (P <0.001) were significantly lower with augmented reality navigation  
 190 system. As far as the diameter of the lesion biopsied, a pathological diagnosis was reached in 96%

191 of cases in group S and 90% of cases in group C. Complete technical success was obtained in all  
192 patients enrolled. Although the diagnostic performance was slightly higher for group S (96%) than  
193 group C (90%), the difference was not statistically significant, probably because PLB is a highly  
194 effective diagnostic procedure that leaves not much room for additional improvement. Finally, the  
195 rate of procedural complications did not show statistically significant differences between subjects  
196 who underwent biopsy assisted by the navigation system with augmented reality (14%) and those  
197 who underwent CT-guided biopsy with the traditional technique (17%). In terms of patients' safety,  
198 40 patients (not included in the sample referred to in the present study) underwent a procedure in  
199 which radiologists reached the target with traditional CT guidance but worked under constant  
200 monitoring of the navigation system with augmented reality. All needle movements and the  
201 trajectory followed by the operator for each individual patient were recorded and compared with the  
202 trajectory that the same operator would have adopted at the beginning of the procedure, following  
203 the instructions provided by the navigation system with augmented reality. In all these cases, and  
204 in subsequent patients included in group S, there were no rough system errors that put patients'  
205 safety at risk or exposed patients to further complications.

206 Grasso and colleagues [Grasso, 2013 (J Int CARS)] performed SIRIO-guided PLB in 197 patients  
207 (121 male, 76 female, mean age  $66.8 \pm 12$  years). Results obtained from this group were compared  
208 with those obtained from a cohort of 72 patients (48 males, 24 females,  $69.1 \pm 10$  years), who  
209 received standard CT-guided PLBs during a 6-month period ahead of SIRIO availability. Differences  
210 in number of CT scans, patient radiation exposure and procedural time in both groups of patients  
211 (SIRIO vs standard CT-guided PLBs) were evaluated by mean and standard deviation. In particular,  
212 number of CT scans was significantly lower in SIRIO group and the 95th percentile was reached  
213 between 8 and 10 scans with SIRIO and between 16 and 18 with the conventional CT approach.  
214 Also, radiation dose was lower in SIRIO group (the 95th percentile was reached between 7.5 and  
215 10mSv with SIRIO; with the standard CT-guided PLBs, the 95th was reached between 20 and  
216 22.5mSv) as well as procedural time. SIRIO-guided PLBs showed a sensitivity of 89.1%, a specificity  
217 of 96.9% and an accuracy of 95.4%. When procedures are stratified in terms of lesion size ( $\leq 20$ mm  
218 and  $> 20$ mm), the accuracy is similar in both groups for lesions  $> 20$ mm (93.7 % in the SIRIO-guided  
219 PLBs group vs. 92.3% in the standard CT-guided PLBs group); on the contrary, SIRIO showed a  
220 marked accuracy improvement for lesions  $\leq 20$ mm (96.8 vs. 91.4%). The overall pneumothorax rate  
221 was significantly lower with SIRIO than in the group undergoing standard CT-guided PLBs (9.1 vs.  
222 16.7%,  $p < 0.01$ ).

223 Iannelli et al [Iannelli, 2018] studied 200 prospective patients from January 2014 to December 2015  
224 to evaluate the effectiveness of SIRIO in performing lung biopsies, with particular attention to lesions  
225 smaller than 1 cm, compared to traditional procedure. Results have been divided according to lesion  
226 dimensions. The authors found a statistically significant decrease in the execution of biopsies with  
227 lesions  $\leq 10$ mm with greater efficacy in terms of reduction of the absorbed dose (mean 43.4 vs  
228 131.5, SD 0.9 vs 5.1), procedural time (mean 28.9 vs 39.5, SD 0.8 vs 4.6) and of the CT scanner  
229 tube (mean 17.5 vs 35.4, SD 1.0 vs 4.2), compared to the standard procedures. Regarding safety  
230 issues, using SIRIO safest and shorter path can be there followed, with the option of avoiding  
231 adjacent anatomical structures. This has sharply reduced peri-procedural complications compared  
232 to standard practice (9 vs 5).

233 The article of Grasso et al [Grasso, 2013 (Eur Radiol)] studied 52 consecutive patients. In a previous  
234 experience [Grasso, 2013 (J Int CARS)], authors validated an optical CT based navigation system  
235 while performing percutaneous lung biopsies (PLBs). The aim of the present study was to evaluate  
236 the optical navigation system while performing PLBs with a low-dose (LD) CT protocol in terms of  
237 the technical feasibility, patients' radiation exposure, imaging quality and complication rate. Patients  
238 were enrolled according to the intention-to-treat principle and randomized into two different groups:  
239 group 1 (15 male, 10 female) underwent PLBs under guidance of the optical CT-based navigation  
240 system; group 2 (16 male, 11 female) underwent PLBs under guidance of the optical CT-based  
241 navigation system with an LD CT protocol. A two-sample t-test was performed to compare the  
242 patients' demographics (age), lesion size (similar in both groups), chest radiation dose (significantly

243 lower in group 2), needle path length (similar in both groups), needle repositioning and procedural  
244 time (similar in both groups);  $P < 0.05$  was considered significant. Descriptive statistics were used for  
245 technical success (100% in both groups), lesion location, diagnosis, complication rate (similar in  
246 both groups) and group 2 CT image quality (in group 2 was always rated as adequate and as  
247 excellent in 15 cases (56.0 %)). The radiation dose to the patients' chest was significantly lower in  
248 group 2. The needle path was similar in both groups as well as needle and procedural time.  
249 Complication rates were similar in both groups.

#### 250 *Non comparative studies*

251 Faiella et al [Faiella, 2018] enrolled a total of 496 patients (mean age  $69.4 \pm 9.6$ ; 307 male, 189  
252 female) from July 2012 to March 2016 in order to investigate and validate the CT navigation system  
253 and to evaluate PLB accuracy based on dimension and location of suspected lesions. Evaluation of  
254 maximum lesion diameter (LD), lesion localization, procedural time (PT), histological sample validity,  
255 minimum lesion distance from pleural surface (DPS), needle distance travelled during the procedure  
256 (DTP) and recovery requirement for major complications were recorded for each patient. The  
257 radiation dose to the patient's chest was estimated by means of the total dose-length product (TDLP)  
258 and then the effective dose was obtained by a specific formula. All patients underwent PLB. In all  
259 cases a CT image demonstrating the needle tip within the target lesion was obtained. In each patient,  
260 the procedure was successfully completed, with a definitive histopathological diagnosis in 96.2% of  
261 cases. Complications were recorded in 156 patients (31.4%): in 133 patients, were minor (26.8%),  
262 in only 23 patients (4.6%) were recorded major complications. The use of the navigation system  
263 allowed a good success in terms of diagnostic value, even in small lesions and in case of wide DTP.  
264 About procedural complications, this study demonstrated that in patients with major complications  
265 (31.4%) there was a higher value of DTP compared to patients with minor complications. Procedural  
266 time (mean 29.5 min) decreased compared to previous published studies adopting standard CT-  
267 guidance (up to 60 min). It has been assessed the difference between a CT-guided lung biopsy  
268 performed with or without the interposition of a bone structure, which can be avoided using the  
269 navigation system: the difference between DPS and DTP was statistically significant ( $p < 0.01$ ) both  
270 in patients presenting an interposed structure and in the whole population. These properties proved  
271 to be more useful in difficult cases, particularly when the pulmonary lesion was small ( $< 15$  mm) or  
272 located at a greater distance from the pleural margin ( $> 45$  mm) or close to critical structures.

273 In their article, Caparelli et al [Caparelli, 2015] retrieved evidence in terms of effectiveness from  
274 other comparative studies [Grasso, 2012; Grasso, 2013 (J Int CARS); Grasso, 2013 (Eur Radiol)],  
275 already cited. As far as safety issues, the first phase of technical evaluation consisted in analyzing  
276 device features that could affect patients' safety. SIRIO H3 did not show any important risks: it has  
277 been designed as aid for minimally invasive surgery, all used materials in the patient area are sterile,  
278 disposable and non-toxic, also not placed in contact with the damaged skin. Following the path  
279 shown on the screen, the operator can insert the needle into the chest and advance it to the target  
280 lesion with a maximum error of 1,5% for the ideal needle trajectory. Dangers associated with the  
281 use of SIRIO are few, due to the fact that is Class I medical device. On the whole, authors evaluate  
282 the risks associated to possible hazard situations as acceptable.

#### 283 **Potential benefits to patients**

284 One of the most interesting aspects of the image-based navigation systems is the use of "augmented  
285 reality", a technique that allows to augment the real visual field of the physician, using information  
286 from a 3-dimensional virtual volume generated with the help of images previously acquired from the  
287 same patient. In this way, radiologists can choose the point of cutaneous insertion more easily and  
288 can rapidly identify the best trajectory to reach the lesion for biopsy, improving also patients' physical  
289 comfort during the procedure [Caparelli, 2015; Grasso, 2012]. The system automatically identifies  
290 the target and graphically traces the trajectory, extending the line on which the needle is localized  
291 in the virtual volume. [Grasso, 2012]. This means the increase of both accuracy (especially for small  
292 lesions) [Hwang, 2018; Grasso 2012; Han, 2018; Faiella, 2018] and safety, as well as the reduction

293 of absorbed ionizing radiation doses [Caparelli, 2015; Grasso, 2012].  
 294 This guiding system has the technical strength of providing real-time virtual navigation of biopsy  
 295 needles based on volume CT data, allowing a flexible selection of the needle route under convenient  
 296 operating conditions. This implies a significant protective factor for pneumothorax because it reduces  
 297 unnecessary redirections [Fior, 2019]. Some highlighted also a reduction in complications rates  
 298 [Faiella, 2018; Iannelli 2018] and a promising procedure in increasing the number of patients and  
 299 nodules to be treated [Caparelli, 2015].  
 300

<input type="checkbox"/> Mortality reduction or increased survival	<input checked="" type="checkbox"/> Reduction of the morbidity	<input type="checkbox"/> Improved quality of life (patient/users)
<input type="checkbox"/> Improved patient monitoring	<input type="checkbox"/> Other: Reduction of the toxicity	<input type="checkbox"/> Not identified

301 **Cost of the technology/procedure**

302 Manufacturers have been contacted through an 'ad hoc' questionnaire sent by e-mail (March 2019).  
 303 As already mentioned, the technology is used on patients who need to undergo a percutaneous CT  
 304 guided procedure, for the diagnosis of a suspected lesion or the therapy of a tumor (ablation type),  
 305 that is neoplastic patients who need to characterize the lesion and/or to treat it. The device is  
 306 intended as assisted technology for biopsy procedures, particularly for pulmonary biopsies. As stated  
 307 by the manufacturer, device cost (VAT excluded) is €139,000. According to manufactures,  
 308 procedures carried out with SIRIO H3 are not reimbursed at national level and services are  
 309 reimbursed according to non-specific DRG, chosen autonomously by each Region (for example, the  
 310 remuneration for percutaneous lung biopsy procedure for each region can range from €250,00 to  
 311 €775,00).

312 Electronic searches to find economic evaluations and cost analysis on SIRIO H3 were performed on  
 313 bibliographic databases (PubMed, and Cochrane Library) in the period April-May 2020. As far as the  
 314 economic aspects, totally the search yielded one paper, that is a preliminary HTA study on guidance  
 315 system for interventional radiology [Caparelli, 2015]. Interventional radiologists using SIRIO H3 into  
 316 clinical practice have been subjected to a questionnaire regarding the use of this technology,  
 317 connected to the patient and also considering their workload. Authors concluded that the costs of  
 318 the technology should be necessarily compared to the provided benefits if the device was covered  
 319 by a public reimbursement policy [Caparelli, 2015]. However, it is worth pointing out that as  
 320 literature is scarce in terms of economic evaluation studies, this lack of evidence is a relevant issue  
 321 that cannot be neglected.  
 322

<input checked="" type="checkbox"/> Increased costs compared to alternative treatments	<input type="checkbox"/> Increased costs due to increased demand	<input type="checkbox"/> Increased costs due to the required investments
<input type="checkbox"/> New costs	<input type="checkbox"/> Other: Reduction of costs linked to the reduction of re-intervention rate	<input type="checkbox"/> Not identified

323

324

325



326 **Potential structural and organisational impact**

327 ***Structural impact***

328 No structural impact is needed for SIRIO H3.

329

<input type="checkbox"/> Increase in requirement of instruments	<input type="checkbox"/> Always be used	<input type="checkbox"/> Can be used only under specific circumstances
<input type="checkbox"/> Decrease in requirement of instruments	<input checked="" type="checkbox"/> Other: no structural impact	<input type="checkbox"/> Not identified

330 ***Organisational impact***

331 Currently, SIRIO H3 is used by an interventional radiologist performing the procedure, a nurse  
 332 handling sterile devices and a radiology technician performing patient CT scan according to a pre-  
 333 set protocol. The manufacturer declared there are no studies investigating learning curve, anyway  
 334 they registered directly from users a considerably reduction of operating times.

335

<input type="checkbox"/> Increase in the number of procedures	<input type="checkbox"/> Re-organisation required	<input checked="" type="checkbox"/> Training required for users
<input type="checkbox"/> Reduction in the number of procedures	<input type="checkbox"/> Other:	<input type="checkbox"/> Not identified

336 **Conclusions**

337 SIRIO H3 is meant for use for all those who have to undergo a biopsy or need drug/ablation therapy  
 338 in a defined area. As reported by the authors, in the retrieved literature, the device appears to be a  
 339 promising alternative to current practice. The available evidence (four comparative studies and two  
 340 non-comparative studies) shows the increase of both accuracy (especially for small lesions) and  
 341 safety, as well as the reduction of absorbed ionizing radiation doses, number of complications and  
 342 procedural time. However, given that no randomized controlled trial is available at present and  
 343 current evidence base is restricted to a few studies, further research is needed to assess the clinical  
 344 properties of the device.

345 **Future prospects**

346 We are not aware of any further potential developments of SIRIO H3 technology. A similar system  
 347 available on market is the CAS-ONE IR (CAScination AG). C-arm cone-beam computed tomography  
 348 (CBCT) systems such as the AXIOM Artis dTA/VB30 (Siemens) and Allura Xper FD20 (Philips  
 349 Healthcare) equipped with virtual guidance (iGuide, Siemens Medical Solutions; XperGuide, Philips  
 350 Healthcare) are possible alternative solutions for percutaneous pulmonary biopsy [Kim, 2015].

351

352 **Evidence searches**

353 Searches of databases including (PubMed, Cochrane Library and ClinicalTrial.gov) in the period  
 354 April-May 2020 using the following keywords:

355

356 **PubMed**

Search		Query	Items found
#6		Search (((((((((((tomography OR x-ray OR spiral))) AND (computed))) OR ((CT OR (CT-guided) OR (computer assisted tomography) OR (Interventional radiology)))))) AND (((((((((((percutaneous OR mininvasiv*)) OR (minimally invasive))) AND ((procedure OR access OR surgery))) AND needle)))))) AND (((Virtual OR (augmented reality)))))) Filters: published in the last 10 years; Humans	43
#5	#3 AND #4	Search (((((((((((((((tomography OR x-ray OR spiral))) AND (computed))) OR ((CT OR (CT-guided) OR (computer assisted tomography) OR (Interventional radiology)))))) AND (((((((((((percutaneous OR mininvasiv*)) OR (minimally invasive))) AND ((procedure OR access OR surgery))) AND needle)))))) AND (((Virtual OR (augmented reality))))))	66
#4		Search ((Virtual OR (augmented reality)))	63874
#3	#2 AND #1	Search (((((((((((((((tomography OR x-ray OR spiral))) AND (computed))) OR ((CT OR (CT-guided) OR (computer assisted tomography) OR (Interventional radiology)))))) AND (((((((((((percutaneous OR mininvasiv*)) OR (minimally invasive))) AND ((procedure OR access OR surgery))) AND needle))))))	3633
#2		Search (((((((((((tomography OR x-ray OR spiral))) AND (computed))) OR ((CT OR (CT-guided) OR (computer assisted tomography) OR (Interventional radiology))))))	839209
#1		Search (((((((((((percutaneous OR mininvasiv*)) OR (minimally invasive))) AND ((procedure OR access OR surgery))) AND needle)	12199

357

358

359 **ClinicalTrials.gov.**

360 No Studies found for: "SIRIO H3" OR "Masmec"

361

362 **Cochrane Library**

363 "virtual navigation system" or "augmented reality" in All Text AND needle in All Text - with  
 364 Cochrane Library publication date Between Jan 2010 and Mar 2020 (Word variations have been  
 365 searched)

366 Items found:

367 *0 Cochrane Reviews matching*

368 *0 Cochrane Protocols matching*

369 **13 Trials matching**

370 *0 Editorials matching*

371 *0 Special collections matching*

372 *0 Clinical Answers matching*

373

374

375

376

377 **Appendix 1**

378 List of included / excluded studies

379

<b>Pubmed, Cochrane, Clinical trials</b>	
Fior D, Vacirca F, Leni D, Pagni F, Ippolito D, Riva L, Sironi S, Corso R. Virtual Guidance of Percutaneous Transthoracic Needle Biopsy with C-Arm Cone-Beam CT: Diagnostic Accuracy, Risk Factors and Effective Radiation Dose. <i>Cardiovasc Intervent Radiol.</i> 2019 May;42(5):712-719. doi: 10.1007/s00270-019-02163-3. Epub 2019 Jan 16	Excluded: intervention does not include SIRIO H3
Iannelli G, Caivano R, Villonio A, Semeraro V, Lucarelli NM, Ganimede MP, Gisone V, Dinardo G, Bruno S, Macarini L, Guglielmi G, Cammarota A. Percutaneous Computed Tomography-Guided Lung Biopsies using a Virtual Navigation Guidance: Our Experience. <i>Cancer Invest.</i> 2018;36(6):349-355. doi: 10.1080/07357907.2018.1498877. Epub 2018 Aug 10.	Included
Hwang EJ, Kim H, Park CM, Yoon SH, Lim HJ, Goo JM. Cone beam computed tomography virtual navigation-guided transthoracic biopsy of small ( $\leq 1$ cm) pulmonary nodules: impact of nodule visibility during real-time fluoroscopy. <i>Br J Radiol.</i> 2018 Jul;91(1087):20170805. doi: 10.1259/bjr.20170805. Epub 2018 Apr 10	Excluded: intervention does not include SIRIO H3
Han Y, Kim HJ, Kong KA, Kim SJ, Lee SH, Ryu YJ, Lee JH, Kim Y, Shim SS, Chang JH. Diagnosis of small pulmonary lesions by transbronchial lung biopsy with radial endobronchial ultrasound and virtual bronchoscopic navigation versus CT-guided transthoracic needle biopsy: A systematic review and meta-analysis. <i>PLoS One.</i> 2018 Jan 22;13(1):e0191590. doi: 10.1371/journal.pone.0191590. eCollection 2018	Excluded: intervention is not virtual navigation
Faiella E, Frauenfelder G, Santucci D, Luppi G, Schena E, Beomonte Zobel B, Grasso RF. Percutaneous low-dose CT-guided lung biopsy with an augmented reality navigation system: validation of the technique on 496 suspected lesions. <i>Clin Imaging.</i> 2018 May - Jun; 49:101-105. doi: 10.1016/j.clinimag.2017.11.013. Epub 2017 Dec 5.	Included
Kim H, Park CM, Lee SM, Goo JM. C-Arm Cone-Beam CT Virtual Navigation-Guided Percutaneous Mediastinal Mass Biopsy: Diagnostic Accuracy and Complications. <i>Eur Radiol.</i> 2015 Dec;25(12):3508-17. doi: 10.1007/s00330-015-3762-8. Epub 2015 Apr 28.	Excluded: intervention does not include SIRIO H3
Mastmeyer A, Hecht T, Fortmeier D, Handels H. Ray-casting based evaluation framework for haptic force feedback during percutaneous transhepatic catheter drainage punctures. <i>Int J Comput Assist Radiol Surg.</i> 2014 May; 9:421-31. doi: 10.1007/s11548-013-0959-7. Epub 2013 Nov 27.	Excluded: not pertinent
Grasso RF, Luppi G, Cazzato RL, Faiella E, D'Agostino F, Beomonte Zobel D, De Lena M. Percutaneous computed tomography-guided lung biopsies: preliminary results using	Included

an augmented reality navigation system. Tumori. 2012 Nov;98:775-82. doi: 10.1700/1217.13503.	
Choo JY, Park CM, Lee NK, Lee SM, Lee HJ, Goo JM. Percutaneous transthoracic needle biopsy of small ( $\leq 1$ cm) lung nodules under C-arm cone-beam CT virtual navigation guidance. . Eur Radiol. 2013 Mar;23:712-9. doi: 10.1007/s00330-012-2644-6. Epub 2012 Sep 14.	Excluded: intervention does not include SIRIO H3
<b>Other sources</b>	
Caparelli C, Carpino G, Brunetti G, Larizza P, Guglielmelli E. A preliminary health technology assessment of a guidance system for interventional radiology. Conf Proc IEEE Eng Med Biol Soc. 2015;2015:450-453. doi:10.1109/EMBC.2015.7318396	Included
Grasso R F, Faiella E, Luppi G, Schena E, Giurazza F, Del Vecchio R, D'Agostino F, Cazzato RL, Zobel BB, "Percutaneous lung biopsy: comparison between an augmented reality CT navigation system and standard CT-guided technique", International Journal of Computer Assisted Radiology and Surgery, 2013, 1-12	Included
Grasso RF, Cazzato RL, Luppi G, et al. Percutaneous lung biopsies: performance of an optical CT-based navigation system with a low-dose protocol. Eur Radiol. 2013;23(11):3071-3076. doi:10.1007/s00330-013-2932-9	Included

380  
381

382  
383

## Bibliography

384

385 ACR–SIR–SPR. PRACTICE PARAMETER FOR THE PERFORMANCE OF IMAGE-GUIDED  
386 PERCUTANEOUS NEEDLE BIOPSY (PNB). 2018. <https://www.acr.org/-/media/ACR/Files/Practice-Parameters/PNB.pdf?la=en> last access 27/06/2019

388 Agenas, Manuale delle procedure HTA,  
389 [http://www.agenas.it/images/agenas/hta/Manuale delle procedure HTA def 20.9.2018 full.pdf](http://www.agenas.it/images/agenas/hta/Manuale_delle_procedure_HTA_def_20.9.2018_full.pdf) ,  
390 last access 24/06/2019

391 Agenas, Piano triennale per la prevenzione della corruzione e la trasparenza,  
392 [http://www.agenas.it/images/agenas/trasparenza/piano\\_corruzione\\_trasparenza/PTPCT\\_2019-2021\\_def.pdf](http://www.agenas.it/images/agenas/trasparenza/piano_corruzione_trasparenza/PTPCT_2019-2021_def.pdf), last access 24/06/2019

394 AIRTUM e AIRTUM Working Group. I numeri del cancro in Italia – 2019. Roma, Settembre 2019.

395 Caparelli C, Carpino G, Brunetti G et al. A preliminary health technology assessment of a guidance  
396 system for interventional radiology. Conf Proc IEEE Eng Med Biol Soc. 2015;2015:450-453.  
397 doi:10.1109/EMBC.2015.7318396

398 Choo JY, Park CM, Lee NK et al. Percutaneous transthoracic needle biopsy of small ( $\leq 1$  cm) lung  
399 nodules under C-arm cone-beam CT virtual navigation guidance. Eur Radiol. 2013 Mar; 23:712-9.  
400 doi: 10.1007/s00330-012-2644-6. Epub 2012 Sep 14.

401 Faiella E, Frauenfelde G, Santucci D et al. Percutaneous low-dose CT-guided lung biopsy with an  
402 augmented reality navigation system: validation of the technique on 496 suspected lesions. Clinical  
403 Imaging 49 (2018) 101–105

404 Fior D, Vacirca F, Leni D et al. Virtual Guidance of Percutaneous Transthoracic Needle Biopsy with  
405 C-Arm Cone-Beam CT: Diagnostic Accuracy, Risk Factors and Effective Radiation Dose. Cardiovasc  
406 Intervent Radiol. 2019 May;42(5):712-719. doi: 10.1007/s00270-019-02163-3. Epub 2019 Jan 16.

407 Grasso RF, Cazzato RL, Luppi G, et al. Percutaneous lung biopsies: performance of an optical CT-  
408 based navigation system with a low-dose protocol. Eur Radiol. 2013;23(11):3071-3076.  
409 doi:10.1007/s00330-013-2932-9

410 Grasso RF, Faiella E, Luppi G et al. Percutaneous lung biopsy: comparison between an augmented  
411 reality CT navigation system and standard CT-guided technique. International Journal of Computer  
412 Assisted Radiology and Surgery, 2013, 1-12

413 Grasso RF, Luppi G, Cazzato RL et al. Percutaneous computed tomography-guided lung biopsies:  
414 preliminary results using an augmented reality navigation system. Tumori. 2012 Nov; 98:775-82.  
415 doi: 10.1700/1217.13503.

416 Han Y, Kim HJ, Kong KA et al. Diagnosis of small pulmonary lesions by transbronchial lung biopsy  
417 with radial endobronchial ultrasound and virtual bronchoscopic navigation versus CT-guided  
418 transthoracic needle biopsy: A systematic review and meta-analysis. PLoS One. 2018 Jan 22;13(1):  
419 e0191590. doi: 10.1371/journal.pone.0191590. eCollection 2018

420 Hwang EJ, Kim H, Park CM et al. Cone beam computed tomography virtual navigation-guided  
421 transthoracic biopsy of small ( $\leq 1$  cm) pulmonary nodules: impact of nodule visibility during real-  
422 time fluoroscopy. Br J Radiol. 2018 Jul;91(1087):20170805. doi: 10.1259/bjr.20170805. Epub 2018  
423 Apr 10

424 Iannelli G, Caivano R, Villonio A et al. Percutaneous Computed Tomography-Guided Lung biopsies  
425 using virtual navigation guidance: our experience. Cancer investigation 2018, vol. 36, no , 349-355.

426 ISTAT. Mortalità per causa - 2017. From I.Stat (<http://dati.istat.it/#> - accessed 09/06/2020).

427 Kim H, Park CM, Lee SM et al. C-Arm Cone-Beam CT Virtual Navigation-Guided Percutaneous

- 428 Mediastinal Mass Biopsy: Diagnostic Accuracy and Complications. *Eur Radiol* (2015) 25:3508–3517.  
429 Mastmeyer A, Hecht T, Fortmeier D, Handels H. Ray-casting based evaluation framework for haptic  
430 force feedback during percutaneous transhepatic catheter drainage punctures. *Int J Comput Assist*  
431 *Radiol Surg*. 2014 May; 9:421-31. doi: 10.1007/s11548-013-0959-7. Epub 2013 Nov 27.  
432

DRAFT

433 **Glossary**

434 **BD/RDM:** Italian Medical device database

435 (<http://www.salute.gov.it/dispositivi/paginainternasf.jsp?id=499&menu=repertorio>).

436 **CND:** Italian medical devices classification (Classificazione Nazionale dei Dispositivi Medici)

437 **CT:** Computed Tomography

438 **FDA:** Food and Drug Administration

439 **MRI:** Magnetic Resonance Imaging

DRAFT