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
**10.4 ANNEX 4 - REPORT FORM FOR FIELD SAFETY CORRECTIVE ACTION**  
**Report Form**  
**Manufacturer's Field Safety Corrective Action Report**  
 Medical Devices Vigilance System  
 (MEDDEV 2.12/1 rev 8)

v.01.13

1. Administrative information	
To which NCA(s) is this report being sent?  Austria, Belgium, Czech Republic, Denmark, Finland, France, Germany, Ireland, Italy, Netherlands, Portugal, Slovenia, Spain, Sweden and United Kingdom.  In addition, this report might be shared with the following competent authorities:  Israel and Kuwait.	
Type of report  <input checked="" type="checkbox"/> Initial report  <input type="checkbox"/> Follow up report  <input type="checkbox"/> Final report	
Date of this report 09-Dec-2016	
Reference number assigned by the manufacturer FCA-2016-062	
FSCA reference number assigned by NCA N/A	
Incidence reference number assigned by NCA	
Name of the co-ordinating national competent authority (if applicable) N/A	
2. Information on submitter of the report	
Status of Submitters  <input checked="" type="checkbox"/> Manufacturer  <input type="checkbox"/> Authorised representative within EEA, Switzerland and Turkey  <input type="checkbox"/> Other (identify the role)	
3. Manufacturer information	
Name Gambro Dasco S.p.A. Medolla Facility	
Contact name Raffaella Turci	
Address Via Modenese 66	
Postcode 41036	City Medolla (MO)
Phone +39.0535.50141	Fax +39 0535 50675
E-mail raffaella_turci@baxter.com	Country Italy
4. Authorized representative information	
Name N/A	
Contact name	
Address	
Postcode	City
Phone	Fax

E-mail	Country
<b>5. National contact point information</b>	
National contact point name Baxter S.p.A	
Name of the contact person Susi Soatti	
Address Via Modenese,66	
Postcode 41036	City Medolla (MO)
Phone +39 053550118	Fax +39 0053550833
E-mail Susi_soatti@baxter.com	Country Italy
<b>6. Medical device information</b>	
Class  <input type="checkbox"/> AIMD Active implants <input type="checkbox"/> MDD Class III <input type="checkbox"/> MDD Class IIb <input checked="" type="checkbox"/> MDD Class IIa <input type="checkbox"/> MDD Class I	<input type="checkbox"/> IVD Annex II List A <input type="checkbox"/> IVD Annex II List B <input type="checkbox"/> IVD Devices for self-testing <input type="checkbox"/> IVD General
Nomenclature system (preferable GMDN) GMDN	Nomenclature code 34999
Nomenclature text : Haemodialysis tubing set, single-use	
Commercial name/brand name/make 1) ArtiSet HD SN HC 2) Artiset HD DNL HC	
Model number  1) 114533 2) 955364	Catalogue number  N/A
Serial number(s)  N/A	lot/batch number(s) 1) Lot: 1634;1635;1636;1637;1638;1639;1640;1641;1642;1644. 2) Lot: 1611;1612;1613;1614;1615;1617;1618;1619;1620;1621;1622; 1623;1624;1625;1626;1627;1628;1629;1630;1631;1634;1635;1636; 1637;1638;1639;1640;1641;1642;1643.
Device Manufacturing date  N/A	Expiry date N/A
Software version number (if applicable) N/A	
Accessories/associated device (if applicable)	
Notified body (NB) ID- number 0086 – BSI	
<b>7. Description of FSCA</b>	
Background information and reason for the FSCA  Baxter is issuing a Safety Alert after having received complaints reporting disconnection between the Artiset bloodline (luer of arterial and/or venous patient connector) and the patient access site (needle/catheter) during treatment. The investigation of disconnection events has determined that the disconnections of the Artiset bloodline from the patient access site were caused by the healthcare provider technique used when connecting the two devices.	
Description and justification of the action (corrective/preventive)  As the disconnection of the bloodline from the patient access site could result in serious adverse health consequences for the patients, such as air embolism and/or blood loss, Baxter is providing its customers with additional instructions on how to ensure a secure connection between the female luer (to patient access) and the male luer (to bloodline), using applicable screwing technique. Besides this, in order to further enforce these recommendations, Baxter is also committed to update the Artiset Bloodlines Instruction For Use, including the additional information for a proper screwing of the patient connectors.	

Advice on actions to be taken by the distributor and the user  
- Distributors are asked to notify their customers of this communication in accordance with their customary procedures  
- Users are asked to:  
1) Follow the additional instructions provided in the Safety Alert communication in order to ensure a proper connection between the Artisect bloodlines and the patient access site.  
2) Acknowledge reception of this Safety Alert (Using the Customer Reply Form)  
3) Forward a copy of this letter to other facilities or departments within their institutions to ensure that those locations are aware of this Important Product Information (where the user is an Hospital/Clinic/Healthcare facility)

Attached please find	FSN Status
<input type="checkbox"/> Field Safety Notice (FSN) in English	<input type="checkbox"/> Draft
<input checked="" type="checkbox"/> FSN in national language	<input checked="" type="checkbox"/> Final
<input type="checkbox"/> Others (please specify):	

Time schedule for the implementation of the different actions  
Implementation date : 14-Dec-2016  
Regional targeted closure date : 06-Oct- 2017  
Local targeted closure date: 22-Sep-2017

These countries within the EEA and Switzerland and Turkey are affected by this FSCA  
Within EEA, Switzerland and Turkey:


<input checked="" type="checkbox"/> AT	<input checked="" type="checkbox"/> BE	<input type="checkbox"/> BG	<input type="checkbox"/> CH	<input type="checkbox"/> CY	<input checked="" type="checkbox"/> CZ	<input checked="" type="checkbox"/> DE	<input checked="" type="checkbox"/> DK	<input type="checkbox"/> EE	<input checked="" type="checkbox"/> ES
<input checked="" type="checkbox"/> FI	<input checked="" type="checkbox"/> FR	<input checked="" type="checkbox"/> GB	<input type="checkbox"/> GR	<input type="checkbox"/> HU	<input checked="" type="checkbox"/> IE	<input type="checkbox"/> IS	<input checked="" type="checkbox"/> IT	<input type="checkbox"/> LI	<input type="checkbox"/> LT
<input type="checkbox"/> LU	<input type="checkbox"/> LV	<input type="checkbox"/> MT	<input checked="" type="checkbox"/> NL	<input type="checkbox"/> NO	<input type="checkbox"/> PL	<input checked="" type="checkbox"/> PT	<input type="checkbox"/> RO	<input checked="" type="checkbox"/> SE	<input checked="" type="checkbox"/> SI
<input type="checkbox"/> SK	<input type="checkbox"/> TR								

Candidate Countries:  
 HR

All EEA, Candidate Countries, Switzerland and Turkey

Others: Israel, Kuwait, Australia, Canada:

8. Comments

I affirm that the information given above is correct to the best of my knowledge.  
  
Signature

Name	City	Date
Susi Soatti	Medolla (MO)	09/12/2016

Submission of this report does not, in itself, represent a conclusion by the manufacturer and/or authorized representative or the national competent authority that the content of this report is complete or accurate, that the medical device(s) listed failed in any manner and/or that the medical device(s) caused or contributed to the alleged death or deterioration in the state of the health of any person

		<b>Busta domestica</b>		<b>myTNT.it</b>	
Con No <b>MY6769 9299</b>		Service		Sender: BAXTER SPA	
Piece 01 of 01 Weight 0,500		Option		P.LE DELL'INDUSTRIA 20 ROMA 00144 RM IT Contact: FABRIZIO TROIANI Tel: 0632491575	
Customer Reference 840061222		Origin Depot RM5		Date 09/12/2016	
Microzona / Fermo Dep.		<b>BT</b>		Delivery Address MINISTERO DELLA SALUTE VIA GIORGIO RIBOTTA 5 00100 ROMA IT Contact: UFFICIO (V) Tel:	
<b>RM21</b>		Dest. Depot <b>RM3</b>		Description of Goods DOCS (C.CAPO) Consignment Volume: 0,001 m <sup>3</sup>	
				<b>BUSTA</b>	
750000223552323372356046					
Special Delivery Instructions OGGETTO: DISPOSITIVO ARTISET - AVVISO DI SICUREZZA					

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