23.11.2018

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For the smallest available stem size (Art. No. 172-916/12) the reaming guide and the guide rod can't be screwed into the intended depth. The reason for this is a reduced depth of the core hole in the stem.

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A reaming guide that is not screwed in completely may lead to a proximal preparation in the depth that is not sufficient and to a neck that can't be positioned correctly.

Corrective reaming can lead to a surgery time extension.

The correct fit of the neck can't be verified by means of the guide rod.

This failure can only occur with the smallest available stem size (Art. No. 172-916/12).

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Safety information and recall of the affected articles from the market

**Misure immediate / Immediate actions:**

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-----

- Should you have any of the Prosthesis stems 172-916/12 in your inventory, please send the products back to Waldemar Link GmbH & Co. KG or confirm implantation.
- An exchange will not take place because this article is taken out of the product portfolio temporarily.
- Affected of this recall are all lot and serial numbers after manufacturing date: 10-2013.
- Please return the fax reply to us in any event until the **12 December 2018** as documentation of the recall. This applies even if you have none of the listed products in stock or if these products do not exhibit the defect in question. Please ensure that all users of the above products within your organization and other relevant persons have been notified of this **safety notice**. If you have transferred the products to third parties, please pass on a copy of this information or notify the contact person indicated below.

The responsible Competent Authorities have been notified of this Field Safety Corrective Action.

**Contatti / Contacts Waldemar Link GmbH & Co. KG:****Gestione reclami / Complaint Management**

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Tel.: +49 40 5 39 95 - 336

[v.pfennig@linkhh.de](mailto:v.pfennig@linkhh.de)

**Responsabile sicurezza / Safety Officer**

Sig. Kai Lemke

Tel.: +49 40 5 39 95 - 142

[vigilance@linkhh.de](mailto:vigilance@linkhh.de)

**Reperto resi / Department of return goods**

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WALDEMAR LINK GmbH & Co. KG