User's Manual

Edi catheter


MAQUET
GETINGE GROUP


TABLE OF CONTENTS

1	Edi catheter	5
2	قسطرة Edi	7
3	Κατετήρ Edi	9
4	Edi-kateter	11
5	Edi kateeter	13
6	سوند ادى	15
7	Edi-katetri	17
8	Cathéter Edi	19
9	Καθετήρας Edi	21
10	Catetere Edi	23
11	Edi カテーテル	25
12	Edi 导管	27
13	Edi 카테터	29
14	„Edi“ kateteris	31
15	Edi-katheter	33
16	Edi-kateter	35
17	Cewnik Edi	37
18	Cateter Aedi	39
19	Cateter Edi	41
20	Κατετερ Edi	43
21	Edi kateter	45
22	Katéter Edi	47
23	Catéter Edi	49
24	Edi-kateter	51
25	สายสวน Edi	53
26	Katetr Edi	55
27	Edi kateteri	57
28	Edi-Katheter	59
29	Edi katéter	61
30	Ống thông Edi	63


T 1

	< 55	< 55	45 - 85	75 - 160	> 140
kg 	0.5 - 1.5	1.0 - 2.0			
Fr/cm	6 / 49	6 / 50	8 / 100	12 / 125	16 / 125 8 / 125

F 1-3B


3a

Fr/cm	Y
16 Fr / 125 cm	$\text{NEX cm} \times 0.9 + 18 = \text{Y cm}$
12 Fr / 125 cm	$\text{NEX cm} \times 0.9 + 15 = \text{Y cm}$
8 Fr / 125 cm	$\text{NEX cm} \times 0.9 + 18 = \text{Y cm}$
8 Fr / 100 cm	$\text{NEX cm} \times 0.9 + 8 = \text{Y cm}$
6 Fr / 50 cm	$\text{NEX cm} \times 0.9 + 3.5 = \text{Y cm}$
6 Fr / 49 cm	$\text{NEX cm} \times 0.9 + 2.5 = \text{Y cm}$

3b

Fr/cm	Y
16 Fr / 125 cm	$\text{NEX cm} \times 0.8 + 18 = \text{Y cm}$
12 Fr / 125 cm	$\text{NEX cm} \times 0.8 + 15 = \text{Y cm}$
8 Fr / 125 cm	$\text{NEX cm} \times 0.8 + 18 = \text{Y cm}$
8 Fr / 100 cm	$\text{NEX cm} \times 0.8 + 8 = \text{Y cm}$
6 Fr / 50 cm	$\text{NEX cm} \times 0.8 + 3.5 = \text{Y cm}$
6 Fr / 49 cm	$\text{NEX cm} \times 0.8 + 2.5 = \text{Y cm}$

1 EDI CATHETER

1.1 INTENDED USE

The Edi catheter is intended for:

- detecting diaphragmatic electrical activity (Edi signals) by which the two ventilator modes NAVA and NIV NAVA can be controlled.
- detecting diaphragmatic electrical activity (Edi signals) by which the respiratory drive from the brain can be monitored.
- administering nutrition, fluids and medications by the naso-gastro-enteric route.

1.2 SAFETY GUIDELINES

WARNINGS!

- The Edi catheter and guide wire are for single use only. Reprocessing will degrade biocompatibility and/or cleanliness. Each Edi catheter may be used for up to 5 days.
- Never re-insert or manipulate the guide wire once the Edi catheter has entered the patient, otherwise the guide wire may penetrate a feeding hole.
- The Edi catheter must be removed from the patient before an MR examination. Due to the metallic electrodes:
 - The strong magnetic and RF fields could generate heat in the electrodes.
 - Artifacts (distortions) may appear in the MR images.

Important:

- The feeding lumen of the Edi catheter is slightly smaller compared to the feeding lumen of a standard feeding tube of the same size; this must be taken into consideration especially for Edi catheters 6 and 8 Fr.
- Do not apply any other substance than water to the Edi catheter. Other substances (lubricants, gels or any other solvents) might destroy the hydrophilic coating and disturb the contact with the electrodes.
- The catheter feeding lumen must be flushed with water after nutrition and medication to prevent clogging.
- Always follow hospital routines for positioning of nasogastric tubes when inserting the Edi catheter.
- The presence of an endotracheal device tends to guide the Edi catheter into the trachea. If any resistance is felt during placement, remove the catheter and reinsert.
- As soon as the Edi catheter position is checked and verified, it is important to attach the Edi catheter securely to the patient.
- If the Edi catheter has been incorrectly positioned, it may detect signals from other muscles than the diaphragm.
- Do not touch the Edi catheter during defibrillation.

Notes:

- Always perform a visual inspection of the equipment before use. Be careful with the Edi catheter electrodes.
- If a guide wire is used, only use a guide wire from MAQUET.
- Continuous movement of the Edi catheter (e.g. against another person's body when holding the baby in kangaroo position) may cause disturbances in the Edi signal.
- External equipment (e.g. heating blankets) may cause disturbances in the Edi signal.
- To prevent contamination of the Edi catheter cable connector, always keep the protective cap on the connector when it is not connected to the ventilator.
- Make sure that the connector cannot hurt the patient; consider distance from patient to connector and possible fixation of connector.
- Follow hospital routines for administration of medication.
- When the feeding lumen is not in use, it is recommended that it is covered in order to avoid incorrect use.
- All disposable parts must be discarded according to hospital routines and in an environmentally safe way.

1.2.1 SYMBOLS

Symbol	Explanation
	Order number
	Number to identify the production batch
	Quantity
	Outer diameter in millimeters
	Inner diameter in millimeters
	Circumference/length
	Use by date
	Do not re-use. Single use only.
	Do not use if packaging is damaged
	Attention – consult documentation
	Keep away from sunlight
	Method of sterilization using irradiation
	Manufacturer
	CE label – indicates compliance with the requirements of the Medical Device Directive 93/42/EEC
	Federal law restricts this device to sale by or on the order of a physician.
	Recycling. Worn-out Edi catheters must be recycled or disposed of properly in accordance with appropriate industrial and environmental standards.
	Phthalate free

1.3 WORKFLOW

- Select the correct Edi catheter size. See table T 1 on page 4.
- Position and calculate insertion distance of Edi catheter. See F 1-3b on page 4.

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For local contact:
Please visit our website
www.maquet.com

0123

GETINGE GROUP is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of ArjoHuntleigh, GETINGE and MAQUET. ArjoHuntleigh focuses on patient mobility and wound management solutions. GETINGE provides solutions for infection control within healthcare and contamination prevention within life sciences. MAQUET specializes in solutions, therapies and products for surgical interventions, interventional cardiology and intensive care.

Edi catheter

User's Manual

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Maquet Critical Care AB