

To the attention of Quality Assurance Dpt or
Regulatory Affairs Dpt or Management

Saint Priest, May 17th 2018

Sujet: **URGENT – FIELD SAFETY NOTICE – RECALL NOTIFICATION LETTER**

Medical device:

ADVANSYS® MLP-DLP; ADVANSYS® TTC; Large QWIX®; TIBIAXYS® and UNI-CP™ - Sterile

Reference:

See list of references in attachment 1

Legal Manufacturer :

*NEWDEAL SAS, Immeuble Séquoïa 2 - 97 allée Alexandre Borodine -
Parc Technologique de la Porte des Alpes - 69800 Saint Priest – France.*

Batch involved:

All non-expired and unused products listed on Appendix 1

Madam, Sir,

Newdeal SAS, a company within Integra LifeSciences Group, has recently identified, through an internal evaluation, the possibility of sealing defect for the plates ADVANSYS® MLP/DLP; ADVANSYS® TTC; TIBIAXYS®; UNI-CP™ and the cortical screws TIBIAXYS® and Large QWIX® packaging. The defect is a non-homogeneous seal or insufficient sealing width and if it were not completely sealed, the sterility of the packaging or the nail itself could be compromised.

Loss of sterility may result in a wound infection that is significant but reversible, requiring intervention beyond standard-of-care. The package defect might not be easily detectable upon visual inspection prior to use. While an adverse health consequence is unlikely to occur based on our health hazard evaluation, Newdeal SAS has made the decision to conduct a voluntary recall of any unused and unexpired products listed on appendix 1.

The review of the available clinical data on the plates ADVANSYS® MLP/DLP; ADVANSYS® TTC; TIBIAXYS®; UNI-CP™ and the cortical screws TIBIAXYS® and Large QWIX® does not raise an abnormal infection rate, consequently no specific follow up for patient implanted is required.

We are notifying you of the recall as our records indicate that you have been supplied with some devices listed on appendix 1.

We kindly ask you to examine your inventory to determine if you have sterile devices listed on appendix 1, please quarantine them.

We also kindly ask you to contact the final customers who may have the affected products and provide them with this letter. If they have affected product, they have to stop using them immediately and remove them from service.

Once the audit of your inventory and your final customers' inventory achieved, please sign and return the "Recall acknowledgment and Return Form" enclosed, by which you confirm that you have received this Recall notification and you intend to fully comply with this Recall notification.

With this form, you will ensure that all the sterile devices ADVANSYS® MLP/DLP; ADVANSYS® TTC; TIBIAXYS®; UNI-CP™ and Large QWIX® affected, will be sent back including those already

Page 1 of 2

Newdeal

Siège Social : Immeuble Séquoïa 2 ■ 97 allée Alexandre Borodine ■ Parc Technologique de la Porte des Alpes ■
69800 Saint Priest ■ France

33 (0)4 37 47 51 51 office ■ 33 (0)4 37 47 51 52 fax ■ integralife.com

Société par Actions Simplifiée au capital de 1.000.000 € ■ NAF 4646Z ■ 412 111 510 RCS Lyon

shipped to your customers. You also confirm that this notification has been forwarded to every concerned customer.

Integra Customer Service will contact you upon receipt of this information to organize the return of the concerned products (Return Merchandise Authorization number assignment).

The receipt of this form ensures that Newdeal has achieved a level of effectiveness in communicating this information.

We also recommend that you keep a copy of this notification and a signed copy of the acknowledgement form for your records.

National Competent Authorities may perform audits of field actions of this nature to verify that our customers have been notified and understand the nature of the field action being taken.

The National Competent Authority of your country has been alerted of this Field Safety Corrective Action.

Please feel free to contact me for any additional questions. Your cooperation is appreciated, and we thank you for your continued collaboration.

Yours Sincerely,



Marilyse Latour
NEWDEAL SAS
Quality Assurance and Regulatory Affairs
Manager

In attached file: Recall acknowledgment and Return Form (2 pages)
Appendix 1 – List of references and batches involved (8 pages)

RECALL ACKNOWLEDGMENT AND RETURN FORM

Medical device:

ADVANSYS® MLP/DLP; ADVANSYS® TTC; Large QWIX®; TIBIAXYS® and UNI-CP™ - Sterile

Reference:

See the list of references in attachment 1

Legal Manufacturer :

**NEWDEAL SAS, Immeuble Séquoia 2 - 97 allée Alexandre Borodine -
Parc Technologique de la Porte des Alpes - 69800 Saint Priest – France.**

Batch involved:

**All non-expired and unused products listed on Attachment 1
May 2018**

Please send the form back to:

By fax/teletype: +33 (0)4 37 47 51 52

Or by e-mail: marilyse.latour@integralife.com

With this form, I confirm that:

I have received, read and understood the information provided in the Integra Recall notification regarding the plates ADVANSYS® MLP/DLP; ADVANSYS® TTC; TIBIAXYS®; UNI-CP™ and the cortical screws TIBIAXYS® and Large QWIX®.

I have transferred this recall letter to the persons to whom I have sold and/or place on consignment the concerned products. I ensure that the form is duly returned to me signed by these persons.

My inventory and my final customers' inventory have been reviewed and the results are as follow (*please tick the appropriate answer*):

Yes, I do have affected product(s) in my inventory or my final customers' inventory. These affected product(s) have been isolated and will be sent back.

Please indicate quantity, lot numbers and circle the reference involved in the table below:

No, I do not have the affected product in my inventory.

I ensure that all the affected products, including those I had already sent to my customers are being quarantined and will be shipped back to Integra.

Distributor / Healthcare facility name

Contact Name

Street Address

City, Country, Postal Code

Telephone

Email

Signature

RECALL ACKNOWLEDGMENT AND RETURN FORM

Medical device:

ADVANSYS® MLP/DLP; ADVANSYS® TTC; Large QWIX®; TIBIAXYS® and UNI-CP™ - Sterile

Reference:

See the list of references in attachment 1

Legal Manufacturer :

**NEWDEAL SAS, Immeuble Séquoia 2 - 97 allée Alexandre Borodine -
Parc Technologique de la Porte des Alpes - 69800 Saint Priest – France.**

Batch involved:

**All non-expired and unused products listed on Attachment 1
May 2018**

Reference	Affected batch	Quantity

APPENDIX 1 – List of references and batches involved by the recall

Reference	Lot	Reference	Lot	Reference	Lot
111740SND	F6H4	111750SND	F70X	111755SND	FBHR
111740SND	F70W	111750SND	F7VP	111755SND	FBRS
111740SND	F9B8	111750SND	F7VQ	111755SND	FCAL
111740SND	F9MY	111750SND	F7VR	111755SND	FEOG
111740SND	FAJ6	111750SND	F9Z0	111755SND	FE0H
111740SND	FATG	111750SND	FA4Y	111755SND	FE0J
111740SND	FBHN	111750SND	FAP1	111755SND	FELL
111740SND	FC9D	111750SND	FB38	111755SND	FELM
111740SND	FD33	111750SND	FBHQ	111755SND	FF4W
111740SND	FD34	111750SND	FC0B	111755SND	FF90
111740SND	FD35	111750SND	FD36	111755SND	FFVS
111740SND	FDQR	111750SND	FD37	111755SND	FFVT
111740SND	FDQS	111750SND	FD38	111755SND	FGCU
111740SND	FELG	111750SND	FE0D	111755SND	FH83
111740SND	FEW0	111750SND	FE0E	111755SND	FHCP
111740SND	FF8X	111750SND	FE0F	111755SND	FHCQ
111740SND	FGA8	111750SND	FELK	111760SND	F6H7
111740SND	FGCS	111750SND	FEW2	111760SND	F6H8
111740SND	FH80	111750SND	FF16	111760SND	F7VV
111740SND	FHKY	111750SND	FF4V	111760SND	F7VW
111745SND	F6H5	111750SND	FF8Y	111760SND	F9Z1
111745SND	F8S3	111750SND	FF8Z	111760SND	FA4Z
111745SND	F9R8	111750SND	FFVR	111760SND	FAP2
111745SND	FA4X	111750SND	FG3M	111760SND	FB31
111745SND	FAJ7	111750SND	FG3N	111760SND	FB3C
111745SND	FB30	111750SND	FH82	111760SND	FBHS
111745SND	FBHP	111750SND	FHCM	111760SND	FBWP
111745SND	FBRR	111750SND	FHCN	111760SND	FD39
111745SND	FBWE	111750SND	FHN7	111760SND	FD3A
111745SND	FCOA	111755SND	F6H6	111760SND	FD3B
111745SND	FE0C	111755SND	F7VS	111760SND	FD3C
111745SND	FELH	111755SND	F7VT	111760SND	FE0K
111745SND	FELJ	111755SND	F7VU	111760SND	FE0L
111745SND	FEW1	111755SND	FAJ8	111760SND	FE0M
111745SND	FGCT	111755SND	FATF	111760SND	FE0N

Reference	Lot	Reference	Lot	Reference	Lot
111760SND	FE0P	111770SND	FB32	111790SND	F712
111760SND	FEW3	111770SND	FB3D	111790SND	F7W0
111760SND	FF4X	111770SND	FBWQ	111790SND	F7W1
111760SND	FF91	111770SND	FE0V	111790SND	F7W2
111760SND	FF92	111770SND	FE0W	111790SND	F9RS
111760SND	FFVU	111770SND	FE0X	111790SND	F9Z6
111760SND	FFVV	111770SND	FELN	111790SND	FA56
111760SND	FFVW	111770SND	FEW4	111790SND	FB2C
111760SND	FGCV	111770SND	FF4Z	111790SND	FB3A
111760SND	FHCR	111770SND	FFVZ	111790SND	FBWK
111765SND	F6H9	111775SND	F6HA	111795SND	F6HC
111765SND	F7VX	111775SND	F70Z	111795SND	F7W3
111765SND	F7VY	111775SND	F9B9	111795SND	F7W4
111765SND	F7VZ	111775SND	F9Z3	111795SND	F9EA
111765SND	F9Z2	111775SND	FA53	111795SND	F9RT
111765SND	FA50	111775SND	FB39	111795SND	F9Z7
111765SND	FAP3	111775SND	FBWG	111795SND	FA57
111765SND	FB29	111775SND	FF50	111795SND	FBWL
111765SND	FBHT	111780SND	F61H	111795SND	FE0Y
111765SND	FBWF	111780SND	F710	111795SND	FE0Z
111765SND	FD3D	111780SND	F8S4	111795SND	FELP
111765SND	FD3E	111780SND	F9RQ	111795SND	FEW5
111765SND	FD3F	111780SND	F9Z4	111795SND	FF52
111765SND	FE0Q	111780SND	FA54	111800SND	F61K
111765SND	FE0R	111780SND	FB2A	111800SND	F7W5
111765SND	FE0S	111780SND	FB33	111800SND	F7W6
111765SND	FE0T	111780SND	FBWH	111800SND	F9R9
111765SND	FE0U	111785SND	F6HB	111800SND	F9RA
111765SND	FF4Y	111785SND	F711	111800SND	F9Z8
111765SND	FFVX	111785SND	F8S5	111800SND	FA58
111765SND	FFVY	111785SND	F9N1	111800SND	FAP4
111765SND	FHCU	111785SND	F9RR	111800SND	FB2D
111765SND	FHCV	111785SND	F9Z5	111800SND	FB3B
111770SND	F6QQ	111785SND	FA55	111800SND	FBWM
111770SND	F70Y	111785SND	FB2B	111800SND	FDQT
111770SND	F9E9	111785SND	FBWJ	111800SND	FDQU
111770SND	FA51	111785SND	FBWJ/1	111800SND	FELQ
111770SND	FA5L	111785SND	FHCW	111800SND	FELR

Reference	Lot	Reference	Lot	Reference	Lot
111800SND	FEW6	111815SND	F9EB	121745SND	FB2S
111800SND	FFW3	111815SND	FAJ9	121745SND	FBJ6
111800SND	FGCW	111815SND	FBEA	121745SND	FC2J
111805SND	F5QC	111815SND	FBRT	121745SND	FD84
111805SND	F6HD	111815SND	FD3K	121745SND	FE1H
111805SND	F8S6	111815SND	FD3L	121745SND	FE1J
111805SND	F9PV	111815SND	FE12	121745SND	FFX8
111805SND	F9Z9	111815SND	FELW	121745SND	FHD9
111805SND	FA59	111815SND	FEW7	121750SND	F6HN
111805SND	FAP5	111815SND	FF57	121750SND	F7WK
111805SND	FBWN	111815SND	FFW6	121750SND	F8YB
111805SND	FD3G	111815SND	FGCY	121750SND	F9U8
111805SND	FE10	111820SND	E62K/G	121750SND	FA4M
111805SND	FELS	111820SND	F7W9	121750SND	FAP7
111805SND	FELT	111820SND	F9EC	121750SND	FB2H
111805SND	FF55	111820SND	F9RB	121750SND	FBJ7
111805SND	FF93	111820SND	FBHV	121750SND	FC2K
111805SND	FFW4	111820SND	FC9E	121750SND	FD85
111805SND	FH4U	111820SND	FD3M	121750SND	FD86
111805SND	FH4V	111820SND	FD3N	121750SND	FE1K
111805SND	FHLO	111820SND	FDQV	121750SND	FE1L
111810SND	E62H/G	111820SND	FELX	121750SND	FFX9
111810SND	F7W7	111820SND	FELY	121750SND	FFXA
111810SND	F9PW	111820SND	FELZ	121750SND	FGD5
111810SND	FATE	111820SND	FGAA	121755SND	F6HP
111810SND	FBHU	121740SND	F7WH	121755SND	F7WL
111810SND	FC0C	121740SND	FAJD	121755SND	FAJF
111810SND	FD3H	121740SND	FB2R	121755SND	FB2J
111810SND	FD3J	121740SND	FBJ5	121755SND	FBJ8
111810SND	FE11	121740SND	FD83	121755SND	FBRV
111810SND	FELU	121740SND	FE1F	121755SND	FC2L
111810SND	FELV	121740SND	FE1G	121755SND	FD87
111810SND	FF56	121740SND	FEX3	121755SND	FE1M
111810SND	FFW5	121740SND	FF9G	121755SND	FE1N
111810SND	FGA9	121740SND	FGD4	121755SND	FF9H
111810SND	FGCX	121745SND	F7WJ	121755SND	FFXB
111815SND	F7RW	121745SND	F9BC	121755SND	FH53
111815SND	F7W8	121745SND	FAJE	121755SND	FHL1

Reference	Lot	Reference	Lot	Reference	Lot
121760SND	F6HQ	121775SND	F9ZD	121790SND	FFXG
121760SND	F7WM	121775SND	FA4P	121790SND	FH86
121760SND	FAH2	121775SND	FBWU	121795SND	F6HW
121760SND	FB2V	121775SND	FDQW	121795SND	F7BH
121760SND	FB88	121775SND	FDQX	121795SND	FBRW
121760SND	FBJ9	121775SND	FDQY	121795SND	FCCC
121760SND	FD89	121775SND	FDQZ	121795SND	FE1X
121760SND	FE1P	121775SND	FE1U	121795SND	FE1Y
121760SND	FE1Q	121775SND	FFXF	121795SND	FF9K
121760SND	FEX4	121780SND	F6HU	121795SND	FGAC
121760SND	FHDB	121780SND	F7JQ	121795SND	FHDE
121760SND	FHDB/1	121780SND	F9UB	121800SND	F6AR
121760SND	FHDC	121780SND	FA4Q	121800SND	FCAM
121765SND	F6HR	121780SND	FAP8	121800SND	FD8F
121765SND	F8BW	121780SND	FBJB	121800SND	FE1Z
121765SND	F9U9	121780SND	FBWV	121800SND	FF9L
121765SND	FB2T	121780SND	FC2M	121800SND	FF9M
121765SND	FBJA	121780SND	FDR0	121800SND	FH55
121765SND	FD8A	121780SND	FDR1	121800SND	FH87
121765SND	FD8B	121780SND	FDR2	121800SND	FH87/1
121765SND	FE1R	121780SND	FDR2/1	121805SND	F6HX
121765SND	FE1S	121780SND	FHDD/1	121805SND	F9FN
121765SND	FFXC	121785SND	F6HV	121805SND	FCAK
121765SND	FH84	121785SND	F7WN	121805SND	FE20
121770SND	F6HS	121785SND	FAJG	121805SND	FF9N
121770SND	F7WW	121785SND	FBJC	121805SND	FGD6
121770SND	F9RD	121785SND	FD8C	121810SND	F7WP
121770SND	FA4N	121785SND	FD8D	121810SND	F9RE
121770SND	FB2W	121785SND	FD8E	121810SND	FD8G
121770SND	FBWT	121785SND	FE1V	121810SND	FE21
121770SND	FE1T	121785SND	FEX5	121810SND	FEX7
121770SND	FFXD	121785SND	FH54	121810SND	FF9P
121770SND	FH85	121785SND	FH54/1	121810SND	FFXH
121770SND	FH85/1	121790SND	F8YD	121815SND	F7WQ
121775SND	F6HT	121790SND	FB2X	121815SND	FBJD
121775SND	F8YC	121790SND	FCAJ	121815SND	FCCD
121775SND	F9PX	121790SND	FE1W	121815SND	FD8H
121775SND	F9UA	121790SND	FEX6	121815SND	FE22

Reference	Lot	Reference	Lot	Reference	Lot
121815SND	FF9Q	150020SND	FGL1	150120SND	FGM9
121815SND	FF9R	150020SND	FGL2	150120SND	FH7W
121815SND	FGAD	150020SND	FH73	150120SND	FHK2
121820SND	F9BD	150020SND	FH7S	150130SND	EB5W/G
121820SND	FBJE	150020SND	FH7T	150130SND	EDJX/G
121820SND	FD8J	150020SND	FJ4J	150130SND	FAUA
121820SND	FD8J/1	150030SND	FGM5	150130SND	FG0J
121820SND	FD8K	150030SND	FHRR	150130SND	FJ4H
121820SND	FE23	150040SND	F7PB	150200SND	F6RR
121820SND	FF9S	150040SND	F9KG	150200SND	F7X0
121820SND	FG3P	150040SND	FAFE	150200SND	F98B
150010SND	F8X6	150040SND	FBTT	150200SND	FAWZ
150010SND	F9K6	150040SND	FC91	150200SND	FD8L
150010SND	F9YE	150040SND	FGL0	150200SND	FD8M
150010SND	F9YH	150040SND	FH7U	150200SND	FE2T
150010SND	FAC1	150110SND	F97D	150200SND	FEX8
150010SND	FAFD	150110SND	F9KH	150240SND	F6RS
150010SND	FBTV	150110SND	FAU9	150240SND	F88U
150010SND	FC90	150110SND	FCN0	150240SND	F9F1
150010SND	FC93	150110SND	FGGW	150240SND	FBBS
150010SND	FD10	150110SND	FGM7	150240SND	FD8N
150010SND	FGL3	150110SND	FH7V	150240SND	FD8P
150010SND	FGL3/1	150110SND	FJ4G	150240SND	FE2U
150020SND	F6R7	150120SND	F77H	150240SND	FE2V
150020SND	F6R8	150120SND	F77J	150240SND	FE2W
150020SND	F6R9	150120SND	F77K	150240SND	FEX9
150020SND	F8X7	150120SND	F9MJ	150240SND	FF1M
150020SND	F8X8	150120SND	F9YF	150240SND	FFPT
150020SND	F9MH	150120SND	F9YG	150240SND	FGNW
150020SND	FAA6	150120SND	FAA7	150240SND	FGNX
150020SND	FAU7	150120SND	FAZZ	150240SND	FGNY
150020SND	FAU8	150120SND	FBTU	150242SND	F6RU
150020SND	FAZW	150120SND	FC92	150242SND	F88V
150020SND	FAZX	150120SND	FCV9	150242SND	F98F
150020SND	FAZY	150120SND	FFGR	150242SND	FAX0
150020SND	FD11	150120SND	FGGY	150242SND	FCB9
150020SND	FD12	150120SND	FGGZ	150242SND	FE2X
150020SND	FFGP	150120SND	FGM8	150242SND	FE2Y

Reference	Lot	Reference	Lot	Reference	Lot
150242SND	FE2Z	150250SND	FGX4	150270SND	FD8S
150242SND	FF1N	150250SND	FHGG	150270SND	FE39
150242SND	FFPU	150255SND	F6YN	150270SND	FF1S
150242SND	FFPV	150255SND	F6YN/1	150270SND	FG7R
150242SND	FG7K	150255SND	F6YN/2	150270SND	FHGJ
150242SND	FGNZ	150255SND	F6YN/3	150275SND	F6S0
150242SND	FGX2	150255SND	F6YN/4	150275SND	F88Z
150246SND	F6RV	150255SND	F6YN/5	150275SND	F98A
150246SND	F88W	150255SND	F6YN/A	150275SND	F9BW
150246SND	F9F2	150260SND	F6RX	150275SND	FBBX
150246SND	FBBT	150260SND	F7X1	150275SND	FCPS
150246SND	FCBA	150260SND	F8XU	150275SND	FE3A
150246SND	FD8Q	150260SND	F9F5	150275SND	FEXC
150246SND	FE30	150260SND	FAX1	150280SND	F6S1
150246SND	FE31	150260SND	FBBV	150280SND	F7X2
150246SND	FE32	150260SND	FCPR	150280SND	F98D
150246SND	FE33	150260SND	FD8R	150280SND	FAPQ
150246SND	FEXA	150260SND	FE36	150280SND	FBBY
150246SND	FF1P	150260SND	FE37	150280SND	FCPT
150246SND	FFPW	150260SND	FF1R	150280SND	FE3B
150246SND	FG7L	150260SND	FFPY	150280SND	FF1T
150246SND	FGP0	150260SND	FG7P	150280SND	FG7S
150246SND	FGX3	150265SND	F6RY	150285SND	F6S2
150246SND	FH2R	150265SND	F88Y	150285SND	F7X3
150250SND	F6RW	150265SND	F989	150285SND	F9F7
150250SND	F88X	150265SND	F9F6	150285SND	FBBZ
150250SND	F9F3	150265SND	F9VA	150285SND	FD8T
150250SND	FBBU	150265SND	FCBB	150285SND	FE3C
150250SND	FCPQ	150265SND	FDSO	150285SND	FF1U
150250SND	FDSB	150265SND	FE38	150285SND	FG7T
150250SND	FDSC	150265SND	FFPZ	150290SND	EB7X/G
150250SND	FE34	150265SND	FG7Q	150290SND	F6S3
150250SND	FE35	150265SND	FH2T	150290SND	F890
150250SND	FEXB	150270SND	F6RZ	150290SND	F98G
150250SND	FF1Q	150270SND	F8LD	150290SND	FBC0
150250SND	FFPX	150270SND	F9BX	150290SND	FCPU
150250SND	FG7M	150270SND	F9VB	150290SND	FDSE
150250SND	FG7N	150270SND	FBBW	150290SND	FE3D

Reference	Lot	Reference	Lot	Reference	Lot
150290SND	FF1V	181013SND	E8NF/G	181041SND	FOAT
150290SND	FH2U	181013SND	E8NF/G1	181041SND	FOAU
150514SND	EA0T/G	181013SND	E8NZ/G	181042SND	E8P8/G
150514SND	F6W9	181013SND	E8PH/G	181042SND	EB1F/G
150514SND	F7Z1	181013SND	F538	181051SND	FCK9
150514SND	F8EB	181013SND	F5CG	181051SND	FG3H
150514SND	FAXQ	181013SND	F7G7	181051SND	FG3J
150514SND	FCGX	181014SND	E8MD/G	181052SND	E8PA/G
150514SND	FEHH	181014SND	E8NG/G	181052SND	EB1H/G
150514SND	FEXD	181014SND	E8P0/G	330021SND	F8ML
150514SND	FG1C	181014SND	F539	330021SND	FEDJ
150516SND	F6WA	181014SND	F539/1	330021SND	FF38
150516SND	F7Z2	181014SND	F5CH	330023SND	EF05/G
150516SND	F8EC	181021SND	E8NH/G	330023SND	EGNK/G
150516SND	FAY4	181021SND	E8P1/G	330023SND	F8MM
150516SND	FCGY	181021SND	EB18/G	330023SND	F9J2
150516SND	FEHJ	181022SND	E8MF/G	330023SND	FAPZ
150516SND	FEXE	181022SND	E8PQ/G	330023SND	FBE8
150516SND	FG1D	181022SND	EBKU/G	330023SND	FBPB
181001SND	E8PB/G	181023SND	E8MG/G	330023SND	FCFQ
181001SND	F5CA	181023SND	E8P3/G	330023SND	FCJP
181002SND	E8M7/G	181023SND	E8PR/G	330023SND	FEDK
181002SND	E8MR/G	181023SND	F8EG	330023SND	FF39
181002SND	E8PC/G	181031SND	E8N2/G	330023SND	FFH7
181003SND	E8MS/G	181031SND	E8NL/G	330025SND	FAE6
181003SND	E8NB/G	181031SND	E8P4/G	330025SND	FBPC
181003SND	E8NV/G	181031SND	E8PS/G	330025SND	FCFR
181003SND	F537	181031SND	EBKW/G	330025SND	FCJQ
181004SND	E8M9/G	181031SND	EBKW/G1	330025SND	FEDL
181004SND	E8NW/G	181031SND	F8EH	330025SND	FFH8
181004SND	E8PE/G	181032SND	E8P5/G	330030SND	EMGF/1
181011SND	E8PF/G	181032SND	E8P5/G1	330030SND	ENHN
181011SND	F5CE	181032SND	E8PT/G	330030SND	EPQM
181012SND	E8MB/G	181032SND	EB1C/G	330030SND	FBE5
181012SND	E8MV/G	181032SND	EBKX/G	330030SND	FBDP
181012SND	E8PG/G1	181032SND	F8EJ	330030SND	FBPE
181013SND	E8MC/G	181033SND	E8PU/G	330030SND	FBPF
181013SND	E8MW/G	181041SND	FOAS	330030SND	FEDM

Reference	Lot	Reference	Lot	Reference	Lot
330030SND	FEXQ	330220SND	FAQ1	330420SND	F9CK
330030SND	FEXR	330220SND	FBP3	330420SND	FA80
330030SND	FFH9	330220SND	FC2Z	330420SND	FB9L
330030SND	FGHY	330220SND	FD69	330420SND	FBPG
330217SND	F7BY	330220SND	FEDQ	330420SND	FEDU
330217SND	F7BZ	330220SND	FEDR	330420SND	FGJD
330217SND	F7C0	330220SND	FGNC	330425SND	F7MZ
330217SND	F7C1	330220SND	FHLM	330425SND	F8P6
330217SND	F7MV	330225SND	F7N0	330425SND	F8XQ
330217SND	F98K	330225SND	F8P5	330425SND	F9K8
330217SND	F9K7	330225SND	F8WS	330425SND	FB4X
330217SND	FAFY	330225SND	FAQ2	330425SND	FBPH
330217SND	FAQ0	330225SND	FBP4	330425SND	FBPJ
330217SND	FC2Y	330225SND	FEDS	330425SND	FF3A
330217SND	FCJR	330225SND	FHLN	330425SND	FFHD
330217SND	FEDN	330230SND	F7MX	330430SND	F7C5
330217SND	FEDP	330230SND	F8WT	330430SND	F9EU
330217SND	FFHA	330230SND	F9SZ	330430SND	F9SY
330217SND	FFHB	330230SND	FAFZ	330430SND	FA81
330220SND	F7C2	330230SND	FBE6	330430SND	FB9M
330220SND	F7C3	330230SND	FEDT	330430SND	FBE7
330220SND	F7C3/1	330230SND	FEDT/1	330430SND	FBPK
330220SND	F7MW	330420SND	F7C4	330430SND	FEDV
330220SND	F8KW	330420SND	F7MY	330430SND	FFHE
330220SND	F8WR	330420SND	F8WU		
330220SND	F9J1	330420SND	F942		