

**AVVISO URGENTE DI SICUREZZA SUL CAMPO**

**Per Siringhe Portex® Arterial Blood Sampling Line Draw**

<b>Dispositivi interessati:</b>	Siringhe Portex® Arterial Blood Sampling Line Draw
<b>Tipo di Azione:</b>	Azione correttiva di sicurezza sul campo – <b>AVVISO DI SICUREZZA</b>
<b>Data:</b>	24 Aprile 2014
<b>All'attenzione di:</b>	Responsabili della sicurezza, Reparti ospedalieri per patologie respiratorie, infermieri, pronto soccorso, responsabili dei materiali, laboratori diagnostici distributori e altri utilizzatori di questi dispositivi.
<b>Dettagli sui dispositivi interessati:</b>	Allegato 1

Smiths Medical a partire dal 22 Novembre 2013 ha notificato ai suoi clienti una Azione Correttiva sul Campo riguardante alcune siringhe Portex® Arterial Blood Sampling Line Draw (siringhe ABL), descritte nell'Allegato 1. Il Ministero della Salute è stato informato di questo Avviso di Sicurezza.

**Nella notifica iniziale NON tutti i lotti coinvolti erano stati identificati. Vi preghiamo di leggere l'Allegato 1 per una lista COMPLETA di tutti i lotti coinvolti e di agire di conseguenza.**

Smiths Medical ha identificato un problema riguardante le siringhe ABL, quando utilizzate con l'analizzatore di gas sanguigni di marca Radiometer®. La punta dello stantuffo della siringa non rimane in posizione quando la sonda dell'analizzatore entra nella siringa. Questo causa delle difficoltà durante l'aspirazione del campione o può causare l'espulsione dello stantuffo dalla siringa. Smiths Medical non ha ricevuto nessuna segnalazione di danno come conseguenza di questo problema.

Gli strumenti Radiometer® per aspirare i campioni di sangue utilizzano una loro unica tecnologia che si basa su l'attrito tra il pistone e la siringa per arrestare la sonda dell'apparecchiatura. Questo problema non è stato identificato in nessun altro analizzatore di gas sanguigni.

**Solo i codici e i lotti elencati nell'Allegato 1 sono coinvolti in questo Avviso di Sicurezza.**

**Azioni che dovrà intraprendere l'utilizzatore:**

Se non utilizzate uno strumento Radiometer®, il prodotto rimane sicuro ed efficace per essere utilizzato secondo la sua destinazione d'uso.

**Utilizzate uno strumento Radiometer®:**

1. Esaminate il vostro inventario e isolate immediatamente le siringhe elencate nell'Allegato 1.
2. Completate e restituite il modulo di conferma (Allegato 2) tramite fax al **02/90165360** o tramite e-mail a [ordini.italia@smiths-medical.com](mailto:ordini.italia@smiths-medical.com), entro cinque giorni dal ricevimento di questo Avviso di Sicurezza. Sarete forniti di istruzioni complete su come restituire il prodotto coinvolto per la sostituzione.

**Non utilizzate uno strumento Radiometer®:**

Completate e restituite il modulo di conferma (Allegato 2) tramite fax al **0773/4084350** o tramite e-mail a [ordini.italia@smiths-medical.com](mailto:ordini.italia@smiths-medical.com), entro cinque giorni dal ricevimento di questo avviso di sicurezza.

Nessun'altra azione è richiesta in quanto il prodotto rimane sicuro ed efficace per essere utilizzato secondo la sua destinazione d'uso.

**Trasmissione di questo Avviso di Sicurezza Urgente:**


Questo avviso deve essere inoltrato a chiunque all'interno della vostra organizzazione possa esserne interessato, compresi gli utilizzatori finali o altre organizzazioni a cui i dispositivi coinvolti potrebbero essere trasferiti. Se avete distribuito i prodotti coinvolti ad altre persone o strutture, siete pregati di inoltrare immediatamente una copia di questo Avviso Urgente di Sicurezza sul Campo.

Siete pregati di mantenere la vostra attenzione su questo Avviso e le sue conseguenze per un periodo di tempo appropriato in modo da assicurare l'efficacia dell'azione correttiva.

Per qualsiasi prodotto relativo a questi prodotti o per qualsiasi domanda, per favore contattate il Customer Service di Smiths Medical Italia allo **02/90165350**.

Smiths Medical è impegnata a fornire prodotti e servizi di qualità ai suoi clienti. Ci scusiamo per ogni inconveniente che questa situazione potrebbe avere causato.

**Marco Savino**  
**Regional Director, Quality Systems**  
**Smiths Medical**  
Tel: +39 0773 4084810  
marco.savino@smiths-medical.com



ALLEGATO 2

**AVVISO URGENTE DI SICUREZZA SUL CAMPO  
MODULO DI CONFERMA**

**Per Siringhe Portex® Arterial Blood Sampling Line Draw**

Si prega di compilare e restituire questo modulo tramite fax al numero: **02/90165360** alla Sig.ra Cecile Franceschi o inviandone una copia elettronica tramite e-mail a: [ordini.italia@smiths-medical.com](mailto:ordini.italia@smiths-medical.com).

<b>Segnare di seguito tutte le caselle d'interesse:</b>
<input type="checkbox"/> – Abbiamo ricevuto e compreso l'Avviso di Sicurezza. Non utilizziamo lo strumento Radiometer®
<input type="checkbox"/> – Abbiamo del prodotto coinvolto nel nostro inventario (per favore compilate la scheda a pagina 5). Siete pregati di contattarmi per fornirmi le istruzioni necessarie alla restituzione.
<input type="checkbox"/> – Non siamo in possesso del prodotto coinvolto – è stato utilizzato/distrutto.
<input type="checkbox"/> – Non abbiamo più prodotto coinvolto. I dispositivi sono stati trasferiti al seguente indirizzo:

Cliente: \_\_\_\_\_

Departimento : \_\_\_\_\_

Firma: \_\_\_\_\_

Data: \_\_\_\_\_

Indirizzo: \_\_\_\_\_

Indirizzo  
Struttura: \_\_\_\_\_  
\_\_\_\_\_

Indirizzo di  
spedizione: \_\_\_\_\_  
\_\_\_\_\_

Numero di  
Telefono:

Email: \_\_\_\_\_



## Attachment 1 - List of Product Reorder and Lot Numbers

Product Reorder Number	Lot Number	Product Name	Expiration Date
4040-2	2472508	ABG-LINEDRAW 3ML LI M	30-Apr-16
4040-2	2476319	ABG-LINEDRAW 3ML LI M	30-Apr-16
4040-2	2485103	ABG-LINEDRAW 3ML LI M	31-May-16
4040-2	2493959	ABG-LINEDRAW 3ML LI M	31-May-16
4040-2	2498064	ABG-LINEDRAW 3ML LI M	30-Jun-16
4040-2	2511207	ABG-LINEDRAW 3ML LI M	30-Jun-16
4040-2	2521581	ABG-LINEDRAW 3ML LI M	30-Jun-16
4040-2	2525780	ABG-LINEDRAW 3ML LI M	31-Jul-16
4040-2	2525781	ABG-LINEDRAW 3ML LI M	31-Jul-16
4040-2	2533316	ABG-LINEDRAW 3ML LI M	31-Aug-16
4040-2	2547272	ABG-LINEDRAW 3ML LI M	31-Aug-16
4040-2	2574207	ABG-LINEDRAW 3ML LI M	30-Sep-16
4040-2	2578556	ABG-LINEDRAW 3ML LI M	30-Sep-16
4042-2	2463876	ABG-LINEDRAW 3ML LI M	30-Apr-16
4042-2	2467728	ABG-LINEDRAW 3ML LI M	30-Apr-16
4042-2	2467730	ABG-LINEDRAW 3ML LI M	30-Apr-16
4042-2	2472511	ABG-LINEDRAW 3ML LI M	31-May-16
4042-2	2476324	ABG-LINEDRAW 3ML LI M	31-May-16
4042-2	2476325	ABG-LINEDRAW 3ML LI M	31-May-16
4042-2	2480609	ABG-LINEDRAW 3ML LI M	31-May-16
4042-2	2485108	ABG-LINEDRAW 3ML LI M	31-May-16
4042-2	2489875	ABG-LINEDRAW 3ML LI M	31-May-16
4042-2	2493965	ABG-LINEDRAW 3ML LI M	31-May-16
4042-2	2503422	ABG-LINEDRAW 3ML LI M	30-Jun-16
4042-2	2503423	ABG-LINEDRAW 3ML LI M	30-Jun-16
4042-2	2506814	ABG-LINEDRAW 3ML LI M	30-Jun-16
4042-2	2511211	ABG-LINEDRAW 3ML LI M	31-Jul-16
4042-2	2513709	ABG-LINEDRAW 3ML LI M	31-Jul-16
4042-2	2521588	ABG-LINEDRAW 3ML LI M	31-Jul-16
4042-2	2529536	ABG-LINEDRAW 3ML LI M	31-Jul-16
4042-2	2533320	ABG-LINEDRAW 3ML LI M	31-Jul-16
4042-2	2536817	ABG-LINEDRAW 3ML LI M	31-Aug-16
4042-2	2547283	ABG-LINEDRAW 3ML LI M	31-Aug-16
4042-2	2551262	ABG-LINEDRAW 3ML LI M	31-Aug-16
4042-2	2551263	ABG-LINEDRAW 3ML LI M	31-Aug-16
4042-2	2560919	ABG-LINEDRAW 3ML LI M	30-Sep-16
4042-2	2560920	ABG-LINEDRAW 3ML LI M	30-Sep-16
4042-2	2569052	ABG-LINEDRAW 3ML LI M	30-Sep-16
4042-2	2569053	ABG-LINEDRAW 3ML LI M	30-Sep-16
4042-2	2574217	ABG-LINEDRAW 3ML LI M	30-Sep-16
4042E	2459794	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	30-Apr-16

## Attachment 1 - List of Product Reorder and Lot Numbers

Product Reorder Number	Lot Number	Product Name	Expiration Date
4042E	2463880	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	30-Apr-16
4042E	2463881	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	30-Apr-16
4042E	2476326	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	31-May-16
4042E	2480611	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	31-May-16
4042E	2493966	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	31-May-16
4042E	2493967	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	31-May-16
4042E	2506812	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	30-Jun-16
4042E	2506813	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	30-Jun-16
4042E	2529537	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	31-Jul-16
4042E	2529539	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	31-Jul-16
4042E	2533321	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	31-Jul-16
4042E	2536818	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	31-Jul-16
4042E	2536820	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	31-Aug-16
4042E	2547286	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	31-Aug-16
4042LH	2459795	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	30-Apr-16
4042LH	2476328	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	31-May-16
4042LH	2503425	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	30-Jun-16
4042LH	2521590	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	31-Jul-16
4042LH	2529540	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	31-Jul-16
4042LH	2533322	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	31-Aug-16
4042LH	2551265	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	31-Aug-16
4042LH	2560930	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	30-Sep-16
4042LH	2569055	LINE-DRAW, DRY HEPARIN, 3CC SYRINGE, LUER LOCK, FILTER PRO	30-Sep-16
4042LH	2574218	3 mL Line Draw, Luer lock syringe (heparin concentration: 7 I.U. per mL) with Filter-Pro® device	30-Sep-16
4043-2	2459798	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	30-Apr-16

## Attachment 1 - List of Product Reorder and Lot Numbers

Product Reorder Number	Lot Number	Product Name	Expiration Date
4043-2	2463882	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	30-Apr-16
4043-2	2467732	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	30-Apr-16
4043-2	2472512	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	31-May-16
4043-2	2472513	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	31-May-16
4043-2	2476329	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	31-May-16
4043-2	2480613	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	31-May-16
4043-2	2480614	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	31-May-16
4043-2	2489880	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	31-May-16
4043-2	2493968	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	31-May-16
4043-2	2498070	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	30-Jun-16
4043-2	2503427	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	30-Jun-16
4043-2	2506816	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	30-Jun-16
4043-2	2511212	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	31-Jul-16
4043-2	2513711	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	31-Jul-16
4043-2	2521591	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	31-Jul-16
4043-2	2521592	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	31-Jul-16
4043-2	2525786	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	31-Jul-16
4043-2	2529541	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	31-Jul-16
4043-2	2533323	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	31-Jul-16
4043-2	2547290	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	31-Aug-16
4043-2	2547291	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	31-Aug-16
4043-2	2551266	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	31-Aug-16
4043-2	2555870	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	30-Sep-16
4043-2	2560932	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	30-Sep-16
4043-2	2569056	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	30-Sep-16

## Attachment 1 - List of Product Reorder and Lot Numbers

Product Reorder Number	Lot Number	Product Name	Expiration Date
4043-2	2569057	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	30-Sep-16
4043-2	2569058	DRY LINEDRAW, ARTERIAL SAMPLING 3ML L/SLIP SYR, FILT-PRO DEVICE	30-Sep-16
4043-2	2574219	LINE DRAW PLUS DRY, 3CC SYRINGE LUER SLIP, FILTER PRO	31-Oct-16
4043E	2459801	LINE DRAW PLUS DRY, 3CC SYRINGE LUER SLIP, FILTER PRO	30-Apr-16
4043E	2467733	LINE DRAW PLUS DRY, 3CC SYRINGE LUER SLIP, FILTER PRO	30-Apr-16
4043E	2467734	LINE DRAW PLUS DRY, 3CC SYRINGE LUER SLIP, FILTER PRO	30-Apr-16
4043E	2476330	LINE DRAW PLUS DRY, 3CC SYRINGE LUER SLIP, FILTER PRO	31-May-16
4043E	2476331	LINE DRAW PLUS DRY, 3CC SYRINGE LUER SLIP, FILTER PRO	31-May-16
4043E	2485112	LINE DRAW PLUS DRY, 3CC SYRINGE LUER SLIP, FILTER PRO	31-May-16
4043E	2485113	LINE DRAW PLUS DRY, 3CC SYRINGE LUER SLIP, FILTER PRO	31-May-16
4043E	2485115	LINE DRAW PLUS DRY, 3CC SYRINGE LUER SLIP, FILTER PRO	31-May-16
4043E	2493970	LINE DRAW PLUS DRY, 3CC SYRINGE LUER SLIP, FILTER PRO	31-May-16
4043E	2493971	LINE DRAW PLUS DRY, 3CC SYRINGE LUER SLIP, FILTER PRO	31-May-16
4043E	2493972	LINE DRAW PLUS DRY, 3CC SYRINGE LUER SLIP, FILTER PRO	31-May-16
4043E	2506817	LINE DRAW PLUS DRY, 3CC SYRINGE LUER SLIP, FILTER PRO	30-Jun-16
4043E	2511213	LINE DRAW PLUS DRY, 3CC SYRINGE LUER SLIP, FILTER PRO	30-Jun-16
4043E	2525787	LINE DRAW PLUS DRY, 3CC SYRINGE LUER SLIP, FILTER PRO	31-Jul-16
4043E	2529542	LINE DRAW PLUS DRY, 3CC SYRINGE LUER SLIP, FILTER PRO	31-Jul-16
4043E	2536822	LINE DRAW PLUS DRY, 3CC SYRINGE LUER SLIP, FILTER PRO	31-Aug-16
4043E	2538273	LINE DRAW PLUS DRY, 3CC SYRINGE LUER SLIP, FILTER PRO	31-Aug-16
4043E	2542763	LINE DRAW PLUS DRY, 3CC SYRINGE LUER SLIP, FILTER PRO	31-Aug-16
4043G-2	2472514	ABG-LINEDRAW 3ML LI M	30-Apr-16
4043G-2	2489881	ABG-LINEDRAW 3ML LI M	31-May-16
4043G-2	2503429	ABG-LINEDRAW 3ML LI M	30-Jun-16
4043G-2	2513713	ABG-LINEDRAW 3ML LI M	30-Jun-16
4043G-2	2525788	ABG-LINEDRAW 3ML LI M	31-Jul-16
4043G-2	2555872	ABG-LINEDRAW 3ML LI M	31-Aug-16



## Attachment 1 - List of Product Reorder and Lot Numbers

<b>Product Reorder Number</b>	<b>Lot Number</b>	<b>Product Name</b>	<b>Expiration Date</b>
G1777	2475627	ABG UMBILICAL LINEDRAW 3ML LI M	30-Apr-16
G1777	2491431	ABG UMBILICAL LINEDRAW 3ML LI M	31-May-16
G2001	2548332	ABG-LINEDRAW 3ML LI M	31-Aug-16

## Attachment 1 - List of Product Reorder and Lot Numbers

Product Reorder Number	Lot Number	Product Name	Expiration Date
4043W-2	2459802	ABG-LINEDRAW 3ML LI M	30-Apr-16
4043W-2	2472515	ABG-LINEDRAW 3ML LI M	30-Apr-16
4043W-2	2485116	ABG-LINEDRAW 3ML LI M	31-May-16
4043W-2	2489882	ABG-LINEDRAW 3ML LI M	31-May-16
4043W-2	2498071	ABG-LINEDRAW 3ML LI M	31-May-16
4043W-2	2503431	ABG-LINEDRAW 3ML LI M	30-Jun-16
4043W-2	2506818	ABG-LINEDRAW 3ML LI M	30-Jun-16
4043W-2	2511214	ABG-LINEDRAW 3ML LI M	30-Jun-16
4043W-2	2521593	ABG-LINEDRAW 3ML LI M	30-Jun-16
4043W-2	2533326	ABG-LINEDRAW 3ML LI M	31-Jul-16
4043W-2	2542764	ABG-LINEDRAW 3ML LI M	31-Aug-16
4043W-2	2547293	ABG-LINEDRAW 3ML LI M	31-Aug-16
4043W-2	2555873	ABG-LINEDRAW 3ML LI M	31-Aug-16
4043W-2	2560933	ABG-LINEDRAW 3ML LI M	30-Sep-16
4043W-2	2569066	ABG-LINEDRAW 3ML LI M	30-Sep-16
4043W-2	2578565	ABG-LINEDRAW 3ML LI M	31-Oct-16
G1469J	2467804	ABG-LINEDRAW 3ML LI M	31-May-16
G1469J	2467805	ABG-LINEDRAW 3ML LI M	30-Apr-16
G1469J	2480661	ABG-LINEDRAW 3ML LI M	31-May-16
G1469J	2503528	ABG-LINEDRAW 3ML LI M	30-Jun-16
G1469J	2503529	ABG-LINEDRAW 3ML LI M	30-Jun-16
G1469J	2521653	ABG-LINEDRAW 3ML LI M	31-Jul-16
G1469J	2529597	ABG-LINEDRAW 3ML LI M	31-Jul-16
G1469J	2529598	ABG-LINEDRAW 3ML LI M	31-Jul-16
G1469J	2551329	ABG-LINEDRAW 3ML LI M	30-Sep-16
G1469J	2572218	ABG-LINEDRAW 3ML LI M	30-Sep-16
G1469J	2572220	ABG-LINEDRAW 3ML LI M	30-Sep-16
G1725	2530465	ABG-LINEDRAW 1ML/3ML LI	31-Jul-16
G1762	2475168	ABG-LINEDRAW 3ML LI M	30-Apr-16
G1762	2479892	ABG-LINEDRAW 3ML LI M	30-Apr-16
G1762	2483968	ABG-LINEDRAW 3ML LI M	31-May-16
G1762	2492596	ABG-LINEDRAW 3ML LI M	31-May-16
G1762	2496746	ABG-LINEDRAW 3ML LI M	31-May-16
G1762	2500736	ABG-LINEDRAW 3ML LI M	31-May-16
G1762	2505477	ABG-LINEDRAW 3ML LI M	30-Jun-16
G1762	2514016	ABG-LINEDRAW 3ML LI M	30-Jun-16
G1762	2524048	ABG-LINEDRAW 3ML LI M	31-Jul-16
G1762	2529356	ABG-LINEDRAW 3ML LI M	31-Jul-16
G1762	2536681	ABG-LINEDRAW 3ML LI M	31-Jul-16
G1762	2541173	ABG-LINEDRAW 3ML LI M	31-Aug-16
G1762	2554034	ABG-LINEDRAW 3ML LI M	31-Aug-16